

PHILIPS

Hand / wrist
stabilization set

ReCare



PTE7030MA



Added stability for users

Help to get stable and correct hand position

This accessory allows patients with weaker hand or wrist to experience upper body conditioning with a firm grip. With reduction in pressure points, neurological patients with little hand strength can enjoy the benefits of ReCare steppers.

More accessibility for neurological patients

- Glove-like fit secures hand to handle

Highlights

Hand wrap



Padded for comfort, the hand wrap secures hands with hook and loop fastener. It fits like a glove and helps to keep hands in proper positioning.

Specifications

Accessory set

set: 1

Adjustment

-: Hook and loop straps

Product material

-: Neoprene foam rubber

