

PHILIPS

Hand/wrist
stabilisation set

ReCare



PTE7030MA



Added stability for users

Helps to achieve stable, correct hand position

This accessory allows patients with impaired grip strength to experience upper body conditioning with a firm grip. With fewer pressure points, neurological patients with decreased grip strength can enjoy the benefits of ReCare steppers.

More accessibility for neurological patients

- Glove-like fit secures hand to handle

Highlights

Hand wrap



Padded for comfort, the hand wrap secures hands with hook and loop fastener. It fits like a glove and helps to keep hands in proper positioning.

Specifications

Accessory set

set: 1

Adjustment

-: Hook and loop straps

Product material

-: Neoprene foam rubber

