

PHILIPS

Grooming kit

4-in-1



QG3020/10



4-in-1 grooming kit

Create your own look

This Philips professional grooming kit QG3020 with head, beard, nose and ear trimmer is the ultimate styling tool to create the right look for any occasion.

Best performance that feels great

- Contour-following comb for speed and comfort
- SteelWave self-sharpening blades for precise results

Create the look you want

- 9 integrated length settings from 1-18 mm
- Beard and moustache comb with 9 length settings
- Full size trimmer for beard and head hair
- Nose and ear trimmer to trim unwanted hairs

Provides optimal power

- Cordless use for maximum freedom
- Charging stand for convenient storage

Highlights

Contour-following comb

The contour-following comb adjusts to every curve to give fast and comfortable results.

Self-sharpening blades



Revolutionary wave-shaped, self-sharpening blades channel and cut hair for a precise and even result.

9 integrated length settings

Select and locks your desired length within a versatile range of possible lengths.

Beard and moustache comb

Trimming your facial hair exactly how you want.

Full size trimmer

Using the full size trimmer without a comb results in a stubble beard look.

Nose and ear trimmer

A must-have addition to your grooming set to keep your unwanted hairs away.

Cordless use only

Cordless use for maximum freedom

Store and charge stand

Charging stand provides convenient storage and ensures the appliance is fully charged and ready to use

Specifications

Create the look you want

Styling tools: 9 position hair clipper comb

Performance

Trimmer size: 32 mm

Power system

Usage: Corded/cordless, Cordless only

Charging time: 10 hours

Running time: 35 minutes

Maintenance

Lubrication: Blades need no oiling

Cleaning: Cleaning brush

Ease of use

Display: Charging indicator

