

PHILIPS

Grooming kit

Multigroom series
3000

Plus

7-in-1



QG3150/30



Superior visibility, maximum precision

The revolutionary MicroBlade technology is designed to give you unprecedented precision to create the exact look you had in mind.

Best performance that feels great

- MicroBlade technology for unprecedented precision
- Sharper than Titanium*: hardened steel blades for long life

Provides maximum ease of use

- Washable attachments

Create the look you want

- 7 attachments for unlimited styling options


Highlights

7 attachment




A range of attachments to enable you to style, maintain and define the hair and/or facial hair style you like.

MicroBlade technology




The revolutionary MicroBlade Technology is designed to give you unprecedented precision to create the exact look you had in mind.

Sharper than Titanium*



The hardened Inox Steel self-sharpening blades are sharper than Titanium*.

Washable attachments



Washable attachments for fast and convenient cleaning.

Specifications

Cutting system

Number of length settings: 9

Self-sharpening stainless steel blades

Cutting element: Stainless steel blades

Ease of use

Secured length settings

Storage: Stand

Maintenance free - No Oil need

Power system

Charging time: 10 hours

Running time: 35 minutes

Usage: Cordless only

Accessories

Cleaning brush

Service

2-year worldwide guarantee

Create the look you want

Styles: Short beard, Chin straps, Goatee, Moustache, Sideburns, Clip your head hair

Attachments

9 position beard and moustache comb

9 position hair clipper comb

Full size trimmer blade: 32 mm

Micro foil shaver: 21 mm

Nose and ear trimmer

Precision trimmer: 21 mm

