

# PHILIPS

Grooming kit

Pro

9-in-1

Waterproof



QG3190



## Superior visibility, maximum precision

The revolutionary MicroBlade technology is designed to give you unprecedented precision to create the exact look you had in mind.

### Best performance that feels great

- MicroBlade technology for unprecedented precision
- Sharper than Titanium\*: hardened steel blades for long life

### Brings maximum ease of use

- 100% waterproof for shower use and easy cleaning

### Create the look you want

- 9+ attachments to create any style you like from head to toe



# Highlights

## 100% waterproof




100% waterproof for easy use and cleaning.

## MicroBlade technology



The revolutionary MicroBlade Technology is designed to give you unprecedented precision to create the exact look you had in mind.

## Sharper than Titanium\*



The hardened Inox Steel self-sharpening blades are sharper than Titanium\*.

## 9+ attachments



Philips offers you a new grooming kit, that allows you to create any style you like from head to toe with superior precision, giving you that sense of pride of a perfect end result

# Specifications

## Cutting system

Number of length settings: 9  
Self-sharpening stainless steel blades  
Cutting element: Stainless steel blades

## Ease of use

Secured length settings  
Storage: Stand  
Maintenance free - No Oil need  
Wet & Dry: 100 % waterproof for easy use and cleaning  
Cleans under the tap

## Power system

Charging time: 10 hours  
Running time: 35 minutes  
Usage: Cordless only

## Accessories

Cleaning brush

## Service

2-year worldwide guarantee

## Create the look you want

Styles: Short beard, Chin straps, Goatee, Moustache, Sideburns, Clip your head hair, Detailed styling, Eyebrow style

## Attachments

4 position eyebrow comb  
9 position beard and moustache comb  
9 position hair clipper comb  
Full size trimmer blade: 32 mm  
Micro foil shaver: 21 mm  
Micro trimmer: 8 mm  
Nose and ear trimmer  
Precision trimmer: 21 mm

