

PHILIPS

Grooming kit

Multigroom series
3000

Pro

9-in-1

Waterproof



QG3190/00



Superior visibility, maximum precision

The revolutionary MicroBlade technology of this Philips grooming kit is designed to give you unprecedented precision to create the exact look you want without unwanted hair.

Best performance that feels great

- MicroBlade technology for unprecedented precision
- Sharper than Titanium*: hardened steel blades for long life

Provides maximum ease of use

- 100% waterproof for shower use and easy cleaning

Create the look you want

- 9+ attachments to create any style you like from head to toe

Highlights

100% waterproof




100% waterproof for easy use and cleaning.

MicroBlade technology



The revolutionary MicroBlade Technology is designed to give you unprecedented precision to create the exact look you had in mind.

Sharper than Titanium*



The hardened Inox Steel self-sharpening blades are sharper than Titanium*.

9+ attachments



Philips offers you a new grooming kit that allows you to create any style you like from head to toe with superior precision, giving you that sense of pride of a perfect end result

Specifications

Cutting system

Number of length settings: 9
Self-sharpening stainless steel blades
Cutting element: Stainless steel blades

Ease of use

Secured length settings
Storage: Stand
Maintenance free - No Oil need
Wet and Dry: 100 % waterproof for easy use and cleaning
Cleans under the tap

Power system

Charging time: 10 hours
Running time: 35 minutes
Usage: Cordless only

Accessories

Cleaning brush

Service

2-year worldwide guarantee

Create the look you want

Styles: Short beard, Chin straps, Goatee, Moustache, Sideburns, Clip your head hair, Detailed styling, Eyebrow style

Attachments

4 position eyebrow comb
9 position beard and moustache comb
9 position hair clipper comb
Full size trimmer blade: 32 mm
Micro foil shaver: 21 mm
Micro trimmer: 8 mm
Nose and ear trimmer
Precision trimmer: 21 mm

