

PHILIPS


Replaceable blades
last up to 4 months

OneBlade

Trim, edge, shave

1 replaceable blade

Fits on all OneBlade handles



QP210/50



Trim, edge and shave any length of hair Designed to cut hair, not skin

The Philips OneBlade is a revolutionary new hybrid styler that can trim, shave and create clean lines and edges, on any length of hair. Forget about using multiple steps and tools. OneBlade does it all.

Comfortable shave

- Unique OneBlade Technology

Easy to use

- Durable OneBlade

Replaceable blades last up to 4 months

QP210/50

Highlights

Unique OneBlade Technology



The Philips OneBlade has a revolutionary new technology designed for facial styling and body grooming. It can trim, edge and shave any length of hair. Its dual protection system – a glide coating combined with rounded tips – makes shaving easier and more comfortable.

Its shaving technology features a fast-moving cutter (200x per second) so it's efficient – even on longer hairs.

Durable OneBlade



The blades are designed for lasting performance. For optimal performance, you will only need to replace each blade every 4 months*. Replacement is easy and hassle free.

Specifications

Trimming and shaving performance

Shaving system: Contour-following technology, Dual protection system

Trimming system: Contour-following technology

Service

Replacement head QP210, QP220: Replace every 4 months*

OneBlade replacement blade

Replacement blades per package: 1

Fits product type: OneBlade (QP25xx), OneBlade (QP26xx), OneBlade Pro (QP65xx), OneBlade Pro (QP66xx)

