

PHILIPS

Replaceable blade


OneBlade

Trim, edge, shave

2 replaceable blades

Fits on all OneBlade handles

Each lasts up to 4 months*



QP220/50



Trim, edge and shave any length of hair Designed to cut hair, not skin

The Philips OneBlade is a revolutionary new hybrid styler that can trim, shave and create clean lines and edges, on any length of hair. Forget about using multiple steps and tools. OneBlade does it all.

Comfortable shave

- Unique OneBlade technology

Easy to use

- Each replaceable blade lasts 4 months

Highlights

Unique OneBlade technology



The Philips OneBlade has a revolutionary technology designed for facial styling and body grooming. It can trim, edge and shave any length of hair. Its dual protection system – a glide coating combined with rounded tips – makes shaving easier and comfortable. Its

shaving technology features a fast-moving cutter (200x per second) so it's efficient – even on longer hairs.

Durable OneBlade



Durable Philips OneBlade that you only need to replace every 4 months* for optimal performance.

Specifications

Trimming & shaving performance

Shaving system: Contour-following technology, Dual protection system

Trimming system: Contour-following technology

OneBlade replacement blade

Replacement blades per package: 2

Fits product type: OneBlade (QP25xx), OneBlade (QP26xx), OneBlade Pro (QP65xx), OneBlade Pro (QP66xx)

