

PHILIPS

WALITA

Blender

Viva Collection

700 W

ProBlend 6

2.4L jar



RI2134/21



Powerful ProBlend 6 blades technology

Triple action of 6 blades, 750 watt motor and jar

Viva collection blender offers a 750W motor, 2.4L Jar and ProBlend 6 blade combined to create perfect end results for your smoothie and food preparation. It makes your tough blending made easy!

Quality guaranteed

- 2 year guarantee

Easy to use

- Over molded knob for easy grip
- Pulse for mixing and blending
- Water level indicator and handle for easy use

Ease of use

- Dishwasher safe except for the main unit and blade assembly

Supreme blending performance

- ProBlend 6 technology for finer blending
- 2.4L Jar
- 5 speeds for soft and hard ingredients


Highlights

2 year guarantee



With 2 years world-wide warranty.

Over molded knob



Over molded knob for easy grip and better speed control.

Pulse for mixing and blending



You can use "Pulse" mode to control blending degree and time conveniently.

Water level indicator



Water level indicator and handle for easy use and pouring.

ProBlend 6 technology



ProBlend 6 technology for effective blending and mixing

Dishwasher safe



All parts are dishwasher safe, except for the main unit and blade assembly.

2.4L Jar

2.4L Jar for everyday nutrients intake

5 Speeds

Blend at 5 different speeds and to different degrees with a choice of setting.

Specifications

Country of origin

Made in: China

General specifications

Number of speed settings: 5

Product features: Dishwasher safe, Integrated cord storage, Variable speed

Service

2-year worldwide guarantee

Technical specifications

Capacity jar: (max) 2,4 L

Power: 1400 W

Voltage: 200 - 230 V

Design

Color: Red

Finishing

Material jar: Tritan

Material of main body: Metal

