

# PHILIPS

## WALITA

### Multicooker

Daily Collection

3D Heating

5 L

980W



RI3237/76



## Your favorite meals at the touch of a button Risottos, cozidos, pastes, soups and desserts

Free up more time using Philips Walita new Multicookers with smart temperature controls that ensures easy cooking of favorite dishes. It comes with 6 automated cooking programs and manual settings for healthy delicious meals every day.

#### The ultimate in convenience

- Swing handle for easy carrying
- Dishwasher friendly & non-stick inner pot
- Detachable steam vent for easy cleaning
- Front top control panel is easy to use

#### Healthy dishes every day

- 3D heating function ensuring evenly heated dishes
- Inspiring recipe book included
- Thick 2mm pot with 5 layers conducts heat evenly
- 11 preset cooking programs and MyRecipe manual function

#### Smart automatic cooking

- Delayed start for up to 24 hours
- Automatic 24hr keep warm
- Automated cooking programs with smart temperature controls

# Highlights

## Dishwasher friendly



Dishwasher friendly & non-stick inner pot

## Swing handle

Convenient and safe to carry the Philips rice cooker away from the kitchen, or to serve rice in dining room.

## Detachable steam vent

Wash the steam vent frequently to remove residual. Simply remove steam vent from the appliance top lid and clean thoroughly.

## 3D heating function



3D heating function ensuring evenly heated dishes.

## Recipe book included

The inspiring recipe book with healthy, tasty recipes and with expert tips is included

## Delayed start



24 hours easy-to-program preset timer ensures meals are ready on time

## Automatic 24hr keep warm



Automatic 24hr keep warm

## Easy to use panel

Front top control panel is easy to use

## Smart temperature control

Automated programs with smart temperature controls that has optimized heating curves for the cooking of each ingredient or dish.

## Thick 2mm 5-layer pot



Thick 2mm pot with 5 layers conducts heat evenly for healthy delicious dishes.

## 11 preset programs



11 preset cooking programs and MyRecipe manual function for delicious variety of healthy dishes.

# Specifications

## Accessories

Included: Measuring cup, Recipe booklet, Scoop, Spatula, Steaming tray/basket

