

Philips  
Baby Food Thermometer



SBCSC201

## Baby's food always just right!

Digital Baby Food Thermometer you can rely on for your child's health and safety. Gives quick accurate measurement so baby food is never too hot.

### **Fast & Reliable measurement**

- Quick and reliable temperature measurement within seconds

### **Quick & easy food preparation**

- Fits every kind of bottle like oval, round and 'wide-necks'

### **Fast & reliable sterilizing**

- Allows boiling or sterilizing without damaging electronics

**PHILIPS**

# Specifications

## Convenience

- Automatic power off
- Ready indicator

## Accessories

- Included accessories: Battery
- User Manual

## Power

- Battery type: LR41

- Number of batteries: 2
- Operating time on battery: 600 hr
- Removable/ replaceable

## Technical specifications

- Humidity: Waterproof up to 1.5m
- Operating temperature range: 32-45C (90-110F) °C
- Accuracy: +/-0.2C
- Measurement time: 10 s

# Highlights

## Measure within seconds

Never too hot, never too cold - the "Just Right" digital thermometer quickly tells you if the milk or food is safe to eat or drink. It takes the guessing out of warming bottles or jars and is great for checking if food is still warm enough after taking a short break in feeding.

## Fits all bottle and jar types

The product is designed in such a way that all kind of sizes and shapes of bottles and jars fit in.

## Detachable electronics

By detaching the sensitive temperature measuring electronics from the pacifier, you can safely sterilize or boil the part your baby comes into contact with.



Issue date 2014-02-06

Version: 3.0.4

12 NC: 9082 100 01767  
EAN: 87 10895 76082 9

© 2014 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)