



SBCSC280

## Healthy food for babies & toddlers

### 3-in-1 Steaming Centre

Healthy & convenient food preparation for your baby and toddler

#### **Baby's food never too hot**

- Overheat protection stops baby's food ever getting too hot
- Uses the latest PTC heating technology for a safe product
- Temperature control system maintains the ideal temperature

#### **Fast & reliable sterilising**

- Destroy all household bacteria and germs harmful for baby
- Sterilizes 6 regular or 4 wide-neck bottles at a time
- No chemicals needed to disinfect baby's bottles & teats
- Hygienic sterilizing function also ideal for small toys

#### **Healthy baby food**

- Preserves vitamins for healthy food

#### **Quick & easy food preparation**

- Defrost your pre-prepared baby food from the freezer
- Baby's food is safely heated up to 37°C within 5 minutes
- Steams baby's vegetables and fruit for healthy nutrition

# PHILIPS

## Technical specifications

### Convenience

- Automatic audio signal routing
- Automatic power off
- **Fits all bottles:** All baby bottles & jar types

### Power

- **Power supply:** 220-240V

## Product highlights

### Overheat protection

To avoid the danger of accidentally burning your baby's mouth, the automatic switch off ensures the food does not get too hot.

### PTC heating element

The controlled heating process incorporating a PTC heating element and two sensors is ideal for consistently heating food and keeping it warm.

### Temperature control system

The latest heating technology incorporates a PTC heating element and two sensors ensures the water in the beaker quickly reaches the ideal temperature and maintains it at a safe level that won't burn your baby's mouth with food or milk that is too hot.

### Disinfects completely

Protect your newborn baby from any harmful household bacteria or germs by disinfecting her bottles, teats, soothers and small plastic toys.

### Holds 4-6 bottles

Convenient use of the sterilizer means having the disinfected bottles ready when you need them. This means fitting all bottles needed for a day's feeding into the sterilizer in one go.

### No use of chemicals

Sterilization is done by using very hot steam.

### Sterilizes small plastic toys

By disinfecting small plastic toys, you can protect your newborn baby from any harmful household bacteria or germs from playing with their toys.

### Preserves all vitamins

When preparing baby's food, keeping all vitamins is very important for your baby's healthy growth. Steaming vegetables or fruit is the optimum way of making sure these healthy vitamins are perfectly maintained.

### De-frosting baby food

Conveniently defrost any baby food which you prepared in advance and put into the freezer. This means you can have healthy baby food ready whenever your baby needs it.

### Perfect temperature in 5 min.

This bottle warmer reliably heats bottles and jars to exactly the right temperature within three minutes to make sure your hungry baby doesn't need to wait long for their food. And no worries about burning baby's mouth either, thanks to the innovative heating technology.

### Steams all veggies & fruits

Once your baby starts eating solid foods, you can prepare healthy steamed vegetables and fruits while preserving the vitamins and nutritional goodness.

### Fits all bottle and jar types

The product is designed in such a way that all kind of sizes and shapes of bottles and jars fit in.



Date of issue 2006-07-07

Version: 3.0

12 NC: 9082 100 02395  
EAN: 87 10895 79306 3

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips Electronics N.V. or their respective owners.  
© 2006 Koninklijke Philips Electronics N.V.  
All Rights reserved.  
[www.philips.com](http://www.philips.com)