

# PHILIPS

Bath and room  
thermometer

SBCSC550



## Measurement of baby's bath & room temperature

### Digital thermometer

This digital thermometer measures the baby's room & bath temperature. Due to its design, it is also attractive as a toy for baby to play with in the water.

#### **Fast & Reliable measurement**

- Quick and reliable temperature measurement within seconds

#### **Totally safe to play with**

- Complies with all relevant toy and safety standards

# Highlights

## Measure within seconds

Never too hot, never too cold - the "Just Right" digital thermometer quickly tells you if the milk or food is safe to eat or drink. It takes the guessing out of warming bottles or jars and is great for checking if food is still warm enough after taking a short break in feeding.

## Complies with toys standard

The products that comply to the toys standard have been thoroughly tested to ensure they meet these norms and are totally safe.

# Specifications

## Power

**Battery type:** CR44

**Number of batteries:** 2

**Operating time on battery:** 22000 hr

**Removable/ replaceable**

## Technical specifications

**Humidity:** Waterproof up to 1.5m

**Operating temperature range:** 10-40C (50-105F) °C

**Accuracy:** +/-0.5C

**Complies with:** EU & USA Toy regulation

