

PHILIPS

Home Facial Device

VisaPure Advanced

Tailored DualMotion Technology

Cleansing, Massage, Fresh eyes

3 heads, pouch, storage palette

2 intensity settings



SC5370/10



Radiant, revitalized and refreshed skin Advanced skin cleansing and skincare

VisaPure Advanced is a home facial device that brings professional expertise to your skincare routine. Your skin becomes radiant, refreshed and revitalised. Enjoy a clean, soft and revitalized skin with a healthy looking glow.

Easy to use, easily fits in your skincare routine

- The heads are easy to clean with warm water and soap
- Perfectly designed for optimal hygiene

Deep gentle cleansing

- Customized DualMotion with Intelligent Head Recognition
- Increases the absorption of your favorite skincare products

Revitalizing Massage

- Better blood circulation to reveal radiant, revitalized skin
- It feels like 750 gentle finger taps per minute
- Developed with Japanese massage experts

Fresh Eyes

- Gentle massage program with 120 nano-vibrations per second
- Refreshes tired eyes in the morning



asimpleswitch.com

Highlights

Easy to clean heads

The heads are easy to clean. Simply clean them in the sink with warm water and soap.

Customized DualMotion



VisaPure Advanced is equipped with Customized DualMotion technology. Each head that comes with VisaPure will have specified levels of rotation and vibration. Both the handle and the new attachments are equipped with an innovative NFC tag that enables the advanced Intelligent Head Recognition. This means that the handle immediately recognizes the specific head once put on. In this way, you can enjoy dedicated programs of Customized DualMotion technology for different skincare benefits.

Increases absorption of creams



Cleansing with VisaPure means that you remove more make up residue, dead skin cells and dullness. Thanks to the deep cleansing effect, your favorite skincare products such as creams, serums and essences are better absorbed by the skin.


Improves blood circulation



Using the Revitalizing Massage head with its Customized DualMotion program, you will increase the blood circulation and bring life

into your skin. This brings out its natural radiance, leaving it looking revitalized and with a healthy glow.

Relaxes the muscles



Thanks to the elegant and functional design of the head with its five small balls, it feels like 750 gentle finger tapping movements per minute. The Revitalizing Massage program lasts 3 minutes, and you can enjoy a facial massage several times a week.

Developed with massage experts



The Revitalizing massage head has been developed with experts in Japanese facial massage. The dedicated program of DualMotion is inspired by worldwide renowned massage techniques. The massage technique called Petrissage is known worldwide to deliver a deep level massage that stimulates blood circulation and relax the muscles, revealing ultimately radiant and revitalized skin.

120 Nano-vibrations per second



The customized DualMotion program for Fresh Eyes delivers 120 nano-vibrations per second, for a gentle, comfortable and refreshing massage around the eyes. It is specially designed to be gentle for the thin and delicate

skin around the eyes. The Fresh eyes program is a short program of only 30 seconds, so you can quickly use it in your daily morning ritual!

Refreshes tired eyes in 30 sec



The Fresh eyes head is made of a special cool material, with a premium ceramic coating, that gives tired eyes a boost of freshness in the morning.

Designed for optimal hygiene



The silky soft bristle material of all our VisaPure brush heads is a specially chosen, delicate nylon. This helps to ensure the brush heads are super easy to clean and keep fresh.



[asimpleswitch.com](https://www.asimpleswitch.com)

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Benefits

Fresh Eyes: Gentle massage program with 120 nano-vibrations per second

Revitalizing Massage: Feels like 750 fingertips tapping your skin per minute

Skin cleansing: 10X more effective than hand-cleansing*

Microcirculation: Improves microcirculation for a more radiant skin

Gentle: as gentle on your skin as manual cleansing*

Exfoliating: Removes more dead skin cells than manual cleansing*

Absorption: Improves the absorption of your skin care product

Application areas

Face and neck: Cheeks, Chest, Chin, Forehead, Neck, Nose

Technical specifications

Two-way rotation

Customized DualMotion

Luxurious storage palette

Store & dry hygienically: Easily store and dry the heads

Ease of use

Waterproof: Can be used in the shower

Battery indicator: Icon indicates battery life

Easy in your skincare routine: Can be used with topicals

Easy to clean heads: Clean with water and soap

LED indicators: Intensity setting, battery low

Items included

Stand: Charging and storing stand

Storage palette

Travel pouch

Brush heads: Normal skin brush head

Power adapter: 100 - 240 V adapter

Heads included: Fresh eyes, Revitalizing Massage

Power

Voltage: 100-240 V

Power system: Rechargeable battery

Charging time: 6 hours

Service

Warranty: 2-year limited warranty

