

Home Facial Device

VisaPure Advanced

Tailored DualMotion Technology

Deep-pore, precision & massage

3 heads, pouch,storage palette

2 intensity settings



SC5371/10



Reveal visibly pure, refined and revitalized skin With precision cleansing for the nose area

Helps to reduce the appearance of enlarged pores after just one treatment.

3-step routine

- Step 1: Deep pore cleansing
- Step 2: Precision cleansing
- Step 3: Revitalizing Massage

Deep gentle cleansing

- Customized DualMotion with Intelligent Head Recognition
- Increases the absorption of your favorite skincare products
- Unique bristle composition for deep pore cleansing

Easy to use, easily fits in your skincare routine


- The heads are easy to clean with warm water and soap
- Perfectly designed for optimal hygiene

Revitalizing Massage

- It feels like 750 gentle finger taps per minute
- Developed with Japanese massage experts

Highlights

Step 1: Deep pore cleansing



The deep pore brush helps to reduce the appearance of pores. Immediate results shown in 9 out of 10 of women.

Step 2: Precision cleansing



The precision brush is uniquely designed to reach and cleanse the hard to reach areas such as around the nose. It also has a customized program, with more intense rotations and a shorter treatment time dedicated to cleanse the nose area for pure and refined skin.

Step 3: Revitalizing Massage



Using the Revitalizing Massage head with its Customized DualMotion programme, you will increase the blood circulation and bring life into your skin. This brings out its natural radiance, leaving it looking revitalized and with a healthy glow.

Customized DualMotion



VisaPure Advanced is equipped with Customized DualMotion technology. Each head that comes with VisaPure will have specified levels of rotation and vibration. Both the handle and the new attachments are equipped with an innovative NFC tag that enables the advanced Intelligent Head Recognition. This means that the handle immediately recognizes the specific head once put on. In this way, you can enjoy dedicated programs of Customized DualMotion technology for different skincare benefits.

Deep Pore Brush Bristles



The deep pore brush comes with a unique combination of bristles of thicker and thinner bristles. The thin bristles reach the pores and cleanse them well, while the thicker bristles swipe away the dirt. The brush consists of 20,000 bristles smaller than the size of your pores, for a superior cleanse. An easy way to a clean and healthy looking skin without impurities.

Increases absorption of creams



Cleansing with VisaPure means that you remove more make up residue , dead skin cells and dullness. Thanks to the deep cleansing

effect, your favorite skincare products such as creams, serums and essences are better absorbed by the skin.

Relaxes the muscles



Thanks to the elegant and functional design of the head with its five small balls, it feels like 750 gentle finger tapping movements per minute. The Revitalizing Massage program lasts 3 minutes, and you can enjoy a facial massage several times a week.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO₂ emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Benefits

Revitalizing Massage: Feels like 750 fingertips tapping your skin per minute
Skin cleansing: Helps to reduce the appearance of pores in 9 out of 10 of women
Gentle: As gentle on your skin as manual cleansing
Exfoliating: Removes more dead skin cells than manual cleansing
Absorption: Improves the absorption of your skin care product

Application areas

Face and neck: Cheeks, Chest, Chin, Fronthead, Neck, Nose

Technical specifications

Two-way rotation
Customized DualMotion

Luxurious storage palette

Store & dry hygienically: Easily store and dry the heads

Ease of use

Waterproof: Can be used in the shower
Battery indicator: Icon indicates battery life
Easy in your skincare routine: Can be used with topicals
Easy to clean heads: Clean with water and soap
LED indicators: Intensity setting, battery low

Items included

Stand: Charging and storing stand
Storage palette
Travel pouch
Brush heads: Deep Pore brush head, Precision brush head
Heads included: Revitalizing Massage
Power adapter: 100 - 240 V adapter

Power

Voltage: 100-240 V
Power system: Rechargeable battery
Charging time: 6 hours

Service

Warranty: 2-year limited warranty

