

PHILIPS

Accessory head for
VisaPure Advanced

Revitalizing Massage

Intelligent Head

SC6060/00



Radiant and revitalized skin

Stimulates blood flow and relaxes the muscles

The Revitalizing Massage Head is developed with experts in Japanese massage. The customized DualMotion programme increases the blood circulation and relaxes the facial muscles, for a revitalized, radiant skin with a healthy looking glow.

Improves blood circulation for radiant skin

- More oxygen and nutrients are available in the skin

Relaxes the muscles for a revitalizing experience

- Reaches the deeper layers of the skin to relax the muscles

Customized programme of DualMotion technology

- Feels like 750 fingertips tapping your skin per minute

Easy to use, easily fits in your skincare routine

- Can be used with and without creams, lotions and oils
- The Head is easy to clean with warm water and soap
- Use 1-3 times a week for a revitalizing experience

Intelligent Head Recognition, tailored DualMotion

- 3 minutes program, 2 intensity settings, two-way rotation



asimpleswitch.com

Highlights

Improves blood circulation


Using the Revitalizing Massage head with its customized DualMotion programme, you will increase your blood circulation, which will give you a radiant skin with a healthy looking glow. Experience a pleasant, revitalizing sensation on your skin.

Relaxes the muscles



The Revitalizing Massage head provides a deep level massage that reaches the deeper layers of the skin. Thereby it relaxes the small facial muscles and gives you a relaxing, comfortable experience.

750 finger taps per minute



The customized DualMotion programme delivers a revitalizing and relaxing experience. It has been developed with experts in Japanese massage and is designed to mimic

the renowned massage techniques called Petrissage and Tapotement. It can do what your hands cannot do. Thanks to the elegant and functional design of the Head with its five small balls, it feels like 750 fingertips tapping your skin per minute, massaging your skin. The Revitalizing Massage programme lasts 3 minutes, and you can enjoy a facial massage several times a week.

Easy in your skincare routine

Easy to integrate in your skincare routine. The Head can be used with your current skincare products and brands, like creams, serums, oils and essences. You can also use it without any cream if you'd prefer.

Easy to clean

The Head is easy to clean. Simply clean it in the sink with warm water and soap.

3 Minute Massage programme

Thanks to the Intelligent Head Recognition of the VisaPure Advanced range, the device recognizes the Head once it is put on, and automatically triggers a customized programme of DualMotion, with dedicated levels of rotation, vibration, duration of the programme and direction of rotation. The Revitalizing Massage programme uses two-way rotation to be able to follow the lymphatic flow in the skin.

Can be used 1-3 times a week

Use the Revitalizing Massage head 1 to 3 times a week for a relaxing and revitalizing experience.

Specifications

Ease of use

Usage only with: VisaPure Advanced
Easy to clean

Intelligent Head Recognition

Intelligent Heads

Technical specifications

Customized DualMotion

Dedicated rotation, vibration

Developed with experts: Experts in Japanese massage

Intensity settings: 2

Program duration: 3 minutes

Two-way rotation

Rotates left and right



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

