

# PHILIPS

## AVENT

Avent BackPack

Black



SCD138/60



## Sleek, comfortable backpack

Lightweight backpack with stylish design

Designed to keep changing essentials on-hand when you are out and about with your baby. Comfortable, lightweight design with padded yet breathable back and shoulder straps. Easy to organize with quick access to front and main compartments.

### Benefits and features:

- Made of lightweight, wipe clean nylon
- Ensures bottles stay warm or cool
- Easy to organize
- Sleek, slim and versatile design
- Keeps changing essentials on-hand

# Highlights

## Ensures bottles stay warm or cool



3M Thinsulate™ pocket keeps two Avent Bottles of pre-boiled hot water warm or cold ready-mixed formula/breast milk cool for several hours.

## Easy to organize

Enjoy quick, easy access to front and main compartments.

## Sleek, slim and versatile design



Perfect for day trips with baby or just when you're out and about.

## Keeps changing essentials on-hand

Large removable baby changing mat and drawstring laundry bag.

# Specifications

## Country of origin

China

## What is included

- Backpack: 1 pcs
- Baby changing mat: 1 pcs
- Laundry bag: 1 pcs
- Removable strap pocket: 1 pcs

## Development stages

Stages: 0 - 6 months, 6 - 12 months

