

PHILIPS

AVENT

Avent urban bag

Brown



SCD148/50



One bag, two looks

Simply swap the front flap to change the design

The Philips Avent Urban bag SCD148/50 has two stylish front flaps to change the look of the bag with the seasons or your wardrobe and has all the features of an Avent Baby Bag.

Convenience on the go

- Ensures bottles stay warm or cool
- Lightweight, wipe-clean microfiber

Easy to organise

- Wide adjustable shoulder strap
- Choice of two zip-on flaps offers two looks in one
- Detachable pocket for personal items

Highlights

Ensures bottles stay warm or cool



For storage of cold expressed breast milk/formula or pre-boiled hot water for up to 4 hours.

Choice of two zip-on flaps offers two looks in one



Choice of two zip-on flaps offers two looks in one

Lightweight, wipe-clean microfiber



The microfibre is lightweight and easy to clean.

Detachable pocket for personal items



Detachable pocket for personal items

Adjustable shoulder strap

The Philips AVANT Urban bag has a wide adjustable shoulder strap for optimum carrying comfort

Specifications

Country of origin

China

What is included

Baby changing mat: 1 pcs
Laundry bag: 1 pcs
Bag: 1 pcs
Personal items pocket: 1 pcs
Insulated bottle carrier: 1 pcs

Development stages

Stages: 6 - 12 months, 0 - 6 months

Design

Color: Comes in various colors

