

PHILIPS

AVENT

Avent TravelBag

Black



SCD149/60



Stylish, elegant weekender bag

Baby travel bag for longer trips

For weekends away, sleepovers with grandparents or as a carry-on for air travel, this roomy Avent travel bag comes complete with all you need.

Comfort on the go

- Full-size, removable baby changing mat
- Lightweight, wipe-clean microfiber

Easy to organise

- Adjustable and detachable shoulder strap
- Insulated bottle carrier
- Separate bag for your personal items

Highlights

Lightweight, wipe-clean microfiber



The microfibre is lightweight and easy to clean.

Full-size, removable baby changing mat

Full-size, removable baby changing mat

Insulated bottle carrier



Insulated bottle carrier to keep 2 Avent Bottles of pre-boiled hot water warm or cold ready mixed formula/breast milk cool for several hours

Adjustable and detachable shoulder strap

Adjustable and detachable shoulder strap

Separate bag for your personal items

Separate bag for your personal items

Specifications

Country of origin

China

What is included

Baby changing mat: 1 pcs

Travel bag: 1 pcs

Insulated bottle carrier: 1 pcs

Baby laundry bag: 1 pcs

Personal items bag: 1 pcs

Development stages

Stages: 6 - 12 months, 0 - 6 months

Design

Color: Comes in various colors

