

PHILIPS

AVENT

Avent Neoprene
ThermaBag

Mixed



SCD150/11



Feeding on-the-go

Designed to hold feeding bottles, VIA Cups and Magic Cups

The stylish ThermaTote holds 2 Philips Avent feeding bottles, 2 Magic Cups or 4 VIA Cups. Its double insulation layer keeps milk cold or water hot for up to four hours. Lightweight, compact and convenient for travel.

Compatible with complete Philips Avent range

- Designed to hold Avent Bottles, VIA and Magic Cups

Convenience on the go

- Wide adjustable shoulder strap
- Ensures bottles stay warm or cool
- Sleek, slim and versatile design

Highlights

Ensures bottles stay warm or cool



For storage of cold expressed breast milk/formula or pre-boiled hot water for up to 4 hours.

Adjustable shoulder strap

The Philips Avent Urban bag has a wide adjustable shoulder strap for optimum carrying comfort

Designed to hold Avent Bottles, VIA and Magic Cups



Designed to hold Avent Bottles, VIA and Magic Cups

Sleek, slim and versatile design



Perfect for day trips with baby or just when you're out and about.

Specifications

Country of origin

China

What is included

ThermaTote: 1 pcs

Convenience

Travel ease: Holds Avent Bottles and VIA, Holds Avent Magic Cups, Keeps warm or cool for 4 hours

Ease of Use: Compact, Lightweight

Development stages

Stages: 6–12 months, 0–6 months

Design

Colour: Comes in various colours

Material

Neoprene

Maximum comfort

Adjustable shoulder strap

