

# PHILIPS

## AVENT

Avent Neoprene  
ThermaBag

Mixed



SCD150/11



## Feeding on-the-go

Designed to hold feeding bottles, VIA Cups and Magic Cups

The stylish ThermaTote holds 2 Philips Avent feeding bottles, 2 Magic Cups or 4 VIA Cups. Its double insulation layer keeps milk cold or water hot for up to four hours. Lightweight, compact and convenient for travel.

### Compatible with complete Philips Avent range

- Designed to hold Avent Bottles, VIA and Magic Cups

### Convenience on the go

- Wide adjustable shoulder strap
- Ensures bottles stay warm or cool
- Sleek, slim and versatile design

# Highlights

## Ensures bottles stay warm or cool



For storage of cold expressed breast milk/formula or pre-boiled hot water for up to 4 hours.

## Adjustable shoulder strap

The Philips Avent Urban bag has a wide adjustable shoulder strap for optimum carrying comfort

## Designed to hold Avent Bottles, VIA and Magic Cups



Designed to hold Avent Bottles, VIA and Magic Cups

## Sleek, slim and versatile design



Perfect for day trips with baby or just when you're out and about.

# Specifications

## Country of origin

Cambodia

## What is included

ThermaTote: 1 pcs

## Convenience

**Travel ease:** Holds Avent Bottles and VIA, Holds Avent Magic Cups, Keeps warm or cool for 4 hours

**Ease of Use:** Compact, Lightweight

## Development stages

Stages: 6–12 months, 0–6 months

## Design

**Colour:** Comes in various colours

## Material

Neoprene

## Maximum comfort

Adjustable shoulder strap

