

PHILIPS

AVENT

CompactBag

Black



SCD151/60



Easy pack – Easy carry

The essential carry bag for short trips

Designed for a quick trip to the shops or to visit a friend for a couple of hours, this essential stylish bag in premium quality, wipe-clean microfibre, accommodates a feed and change for baby and has an internal pocket for personal items.

Benefits and features:

- Adjustable shoulder strap
- Full-size, removable diaper changing mat
- Main compartment is large enough for clothes, pacifiers
- Made of premium quality, wipe clean microfibre
- Ensures bottles stay warm or cool
- Zippered, lined pocket for keeping wipes moist

Highlights

Ensures bottles stay warm or cool



3M Thinsulate™ pocket keeps hot contents warm, or cold contents cool for several hours.

Specifications

Country of origin

China

What is included

Baby changing mat: 1 pcs

CompactBag: 1 pcs

