

PHILIPS

AVENT

Natural Glass Baby
Bottle Gift Set

Natural

SCD201/01



The most natural way to bottle feed

Natural latch on

This newborn starter giftset includes five Philips Avent Natural Glass bottles, a glass bottle sleeve and two Soothie pacifiers. Our Natural bottle makes it easy to combine breast and bottle feeding.

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple
- Soft nipple designed to mimic the feel of the breast
- Flexible spiral design combined with comfort petals

Pure glass

- Heat resistance
- Pharmaceutical grade glass

Vents air away from your baby's tummy

- Unique anti-colic valve technology

Other benefits

- Bottle sleeve protects 4oz/120ml glass bottles
- One piece, silicone Soothie pacifier helps calm babies

Highlights

Natural latch on



The wide breast shaped nipple promotes a natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Ultra soft nipple



The nipple has an ultra soft texture, designed to mimic the feel of the breast.

Flexible spiral design



Flexible spiral design, combined with our unique comfort petals to create a flexible nipple, allowing for a more natural feed without nipple collapse.

Heat resistance



The Natural glass bottle is heat and thermal shock resistant. It can therefore be safely stored in the fridge, warmed, and is also suitable for sterilization.

Pharmaceutical grade glass

Premium quality borosilicate glass to ensure the highest quality for ultimate purity.

Anti-colic valve technology



Our unique anti-colic valve technology was designed to reduce colic and discomfort, by venting air away from your baby's tummy.

Protective bottle sleeve

The insulating bottle sleeve, allows you to keep the feed at the right temperature and to protect your glass bottle at the same time.

Silicone soothie pacifier



The silicone soothie helps calm your baby. It is used by medical professionals to calm newborns.

Specifications

Material

Bottle: Glass

BPA free*: Bottle made of Borosilicate type 1 glass; Nipple made of silicone; Screw ring made of PP (Polypropylene)

Nipple: Silicone, BPA free*

Sleeve: Oeko-Tex® Class 1 certified, Soft polyester

What is included

8oz Baby bottle: 2 pcs

4oz Baby bottles: 3 pcs

Soothe pacifier: 2 pcs

Glass bottle sleeve: 1 pcs

Ease of use

Bottle use: Easy to assemble, Easy to clean, Easy to hold, Dishwasher & microwave safe

Design

Bottle design: Wide neck, Ergonomic shape

Development stages

Stage: 0 - 6 months

Easy to use

Sleeve use: Easy scale reading, Fits glass bottles, Shock resistant, Warm and cool insulating

Functions

Latch on: Natural latch on, Easy combine breast and bottle

Nipple: Unique comfort petals, Ultra soft and flexible nipple

Anti-colic valve: Advanced anti-colic system

