

PHILIPS

AVENT

Natural All in One
Gift Set

Natural

SCD207/01



The most natural way to bottle feed

Natural latch on

This all-in-one giftset includes all you need to start bottle feeding in a natural way - a bottle warmer, a microwave sterilizer, three Philips Avent Natural bottles, two pacifiers and a bottle brush.

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple
- Soft nipple designed to mimic the feel of the breast
- Flexible spiral design combined with comfort petals

Vents air away from your baby's tummy

- Unique anti-colic valve technology

Other benefits

- Compatible with the Philips Avent Natural range

Other essential accessories included in set

- Microwave sterilizer
- Bottle and baby food warmer
- Orthodontic, symmetrical and collapsible pacifier nipple
- Bottle brush with curved brush head for easy cleaning

Highlights

Natural latch on



The wide breast shaped nipple promotes a natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Ultra soft nipple



The nipple has an ultra soft texture, designed to mimic the feel of the breast.

Flexible spiral design



Flexible spiral design, combined with our unique comfort petals to create a flexible nipple, allowing for a more natural feed without nipple collapse.

Anti-colic valve technology



Our unique anti-colic valve technology was designed to reduce colic and discomfort, by venting air away from your baby's tummy.

Compatibility across the range



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

Microwave sterilizer

Ultra fast and easy to use. Sterilizes 4 bottles in 2 minutes. Eliminates harmful bacteria. If lid is unopened, contents remain sterile for up to 24 hours. Also convenient for travel. 2 min. at 1100-1850 W, 4 min. at 850-1000 W, 6 min. at 500-800 W

Bottle and baby food warmer



Fast and easy to use. Heats gently and evenly. No hot spots so safe for baby.

Orthodontic pacifier

Soft, long-lasting silicone nipple. Orthodontic, flat, symmetrical nipple respects the natural development of palate, teeth and gums.

Curved brush head



Specially curved brush head and molded handle-tip reach the corners of all wide neck bottles, nipples and feeding products for thorough cleaning.

Specifications

Material

Bottle: BPA free*

Nipple: Silicone, BPA free*

What is included

9oz Baby bottle: 1 pcs

4oz Baby bottles: 2 pcs

Microwave sterilizer: 1 pcs

Bottle and baby food warmer: 1 pcs

Pacifier: 2 pcs

Bottle and nipple brush: 1 pcs

Design

Bottle design: Wide neck, Ergonomic shape

Ease of use

Bottle use: Easy to assemble, Easy to clean, Easy to hold, Dishwasher & microwave safe

Development stages

Stage: 0 - 6 months

Functions

Latch on: Natural latch on, Easy combine breast and bottle

Nipple: Unique comfort petals, Ultra soft and flexible nipple

Anti-colic valve: Advanced anti-colic system

