

PHILIPS

AVENT

Classic+ PA
Newborn Starter Set

Newborn Starter Set



SCD210/00



Clinically proven to reduce colic and discomfort*

Easy to clean for perfect hygiene

A handy collection including 4 Classic+ PA bottles (2x 4oz, 2x 9oz), a bottle and teat brush, and a soother 0-6 months. The anti-colic system is integrated into the teat, and the bottles are extra easy to clean and assemble.

Clinically proven to reduce colic and discomfort*

- Anti-colic system proven to reduce colic*
- Fussing is reduced, especially at night

Easy to use, clean and assemble

- Few parts for quick and easy assembly and cleaning
- Leak-free for an enjoyable feeding

High-quality and high-durability


- PA material is light, durable and transparent as glass

Other benefits

- BPA free*
- Orthodontic, symmetrical and collapsible soother teat
- Bottle brush with curved brush head for easy cleaning
- Full compatibility with Philips Avent Classic+ bottle range
- Different teat flow rates available in Classic+ range


Highlights

Less colic



Unlike other bottles, the clinically proven anti-colic system is integrated into the teat, making it easier to assemble the bottle correctly. As your baby feeds, the unique valve on the teat opens to allow air into the bottle instead of your baby's tummy.

Less fussing



Sleep and nutrition are vital to your baby's health and happiness. A randomised clinical trial was carried out to see whether infant baby bottle design affects "infant behavior". The Philips Avent Classic baby bottle was shown to significantly reduce fussing by approximately 28 minutes a day as compared to the comparator bottle (46 min vs 74 min, $p=0.05$). This was especially true during the night-time.*

Easy to assemble and clean



With only 4 parts, a wide bottle neck, and rounded corners, our bottle is easy to quickly assemble and clean thoroughly. Have peace of mind that your bottle is truly clean in no time.

Leak-free

The Philips Avent Classic+ bottle is designed to prevent leakage whilst feeding, for a truly enjoyable feeding experience.

Made of Polyamide

The imported PA material is a premium plastic, light, durable and as transparent as glass.

BPA free



The Philips Avent Classic+ PA bottle is made of BPA free* material (polyamide).

Orthodontic Soother



Flat, drop-shaped symmetrical teats respect the natural development of your baby's palate, teeth and gums, even if the soother ends up upside-down in the mouth.

Curved brush head



Specially curved brush head and molded handle-tip reach the corners of all wide neck bottles, nipples and feeding products for thorough cleaning.

Compatibility across range



This Philips Avent Classic+ PA bottle is compatible with most of the Philips Avent range. We recommend to use the Classic+ bottle with Classic+ teats only

Fits your baby's growing needs



The Philips Avent Classic+ range offers different teat flow rates to keep up with your baby's growth. Remember that age indications are approximate as babies develop at different rates.

Specifications

Country of origin

England

Brush: Made in Malaysia

Bottle

Material: BPA free, Highly durable

Capacity: 2x 4oz/125ml, 2x 9oz/260ml

What is included

Baby bottle: 4 pcs

Bottle and nipple brush: 1 pcs

Soother: 1 pcs

Design

Bottle design: Easy to assemble, Easy to clean, Easy to hold, Wide neck

Functions

Anti-colic: Proven anti-colic system

Development stages

Stages: 0-12 months

