

PHILIPS

AVENT

Newborn Starter Set

SCD234/07



For healthy, active feeding

Clinically proven to reduce colic*

A clinical trial demonstrated that at 2 weeks of age, babies fed with the AVENT Bottle experienced less colic than babies fed with a conventional bottle.


Benefits and features:

- Clinically proven to reduce colic
- Built-in Airflex Valve
- Easy to combine breast and bottle feeding

Highlights

Clinically proven to reduce colic

Duration of reported colic per 24 hours at 2 weeks



A clinical trial demonstrated that at 2 weeks of age, babies fed with the AVENT Bottle experienced less colic than babies fed with a conventional bottle. (www.philips.com/AVENT).

Built-in Airflex Valve



The AVENT Airflex Feeding Bottle uses an Airflex Valve which works with baby's natural feeding action.

Specifications

Country of origin

England

What is included

Airflex Feeding Bottle (260 ml/ 9 oz): 2 pcs
 Airflex Feeding Bottle (125 ml/ 4 oz): 2 pcs
 Extra soft, Slow Flow Teat: 2 pcs
 Extra soft, Newborn Flow Teat: 4 pcs
 Newborn Pacifier: 1 pcs
 Bottle and teat brush: 1 pcs

Development stages

Stages: 0 - 6 months

