

PHILIPS

AVENT

Bottle Feeding
Solutions Set

SCD249/00



Everything you need to sterilise and bottle feed

A convenient feeding set

Includes Electric Steam Steriliser, Electric Bottle Warmer, Thermabag, Feeding Bottles, Teats, Milk Powder Dispenser and more

BPA-Free

- Made from Polyethersulfone (PES)

Electric Steam Steriliser

- Large capacity
- Sterile contents for up to 6 hours

Electric Bottle Warmer

- Heats safely and evenly


Feeding Bottle

- Clinically proven to reduce colic*

Highlights

Clinically proven to reduce colic*

Duration of reported colic per 24 hours at 2 weeks



*A clinical trial demonstrated that at 2 weeks of age, babies fed with the AVENT Bottle experienced less colic than babies fed with a conventional bottle. (www.philips.com/AVENT).

Large capacity



Holds up to six 260 ml/9 oz AVENT Bottles or two Philips AVENT Breast Pumps and accessories.

This product is BPA-Free



BPA is Bisphenol-A. This product uses Polyethersulfone (PES), which is BPA-Free and has a natural honey-coloured tint.

Country of origin

England

What is included

- Electric Steam Steriliser: 1 pcs
- Electric Bottle Warmer: 1 pcs
- ThermaBag: 1 pcs
- Classic Bottle 4 oz: 7 pcs
- Classic Bottle 9 oz: 7 pcs
- Soft, Medium Flow Teat: 6 pcs
- Soft, Variable Flow Teat: 6 pcs
- Soft, Fast Flow teat: 6 pcs
- Milk powder dispenser: 1 pcs
- Magic Cup (200 ml/7 oz): 1 pcs
- Soft spout: 2 pcs
- Magic Handles: 2 pcs
- Cup (180 ml/6 oz): 1 pcs
- Silicone soother: 2 pcs
- Bottle and teat brush: 1 pcs

Development stages

Stages: 0–6 months

