

# PHILIPS

## AVENT

Breastfeeding  
Solutions Set

SCD281/00



## Breastfeeding made easy

### Avent Breast care

Express comfortably, store, feed, and protect your nipples with the Philips Avent Breastfeeding solutions set.

#### **Express**

- Manual breast pump with massage cushion

#### **Store**

- Sealing discs

#### **Feed**

- Bottles clinically proven to reduce colic

#### **Care**

- Moisturizing nipple cream

# Highlights

## Manual breast pump



Being stressed or in a hurry can make it more difficult to release your milk and can also affect your milk supply. So feeling comfortable and relaxed really helps. Our manual breast pump has been designed to increase your comfort when you are pumping.

## Bottles clinically proven



Unique anti-colic system consists of the Philips AVENT nipple and adaptor ring allowing air to flow into the bottle instead of your baby's tummy for better digestion.\*

## Moisturizing nipple cream



Hypoallergenic for mother and baby. Moisturizes and softens dry or sensitive nipples. Easy to apply and no need to remove before breastfeeding.

## Sealing discs



Convert all Philips AVENT bottles into milk storage containers.

# Specifications

## Country of origin

EU

## What is included

- Manual breast pump: 1 pcs
- Classic Bottle 4oz: 3 pcs
- Classic Bottle 9oz: 2 pcs
- Extra soft Newborn Flow Nipple: 2 pcs
- Extra soft Slow Flow Nipple: 2 pcs
- Newborn nipple travel pack: 1 pcs
- Sealing disc for milk storage: 6 pcs
- Moisturizing Nipple Cream: 1 pcs
- Disposable breast pad: 40 pcs
- Bottle and nipple brush: 1 pcs
- Guide to breastfeeding DVD: 1 pcs

## Development stages

Stages: 0 - 6 months

