

PHILIPS

AVENT

Breastfeeding
essentials

Classic

SCD282/00



Breastfeeding made easy

Avent Breast care

Express comfortably, store, feed and protect your nipples with the Philips Avent Breastfeeding essentials set.

Express

- Electronic breast pump remembers your personal rhythm

Store

- Milk storage cups
- Sealing discs

Feed

- Bottles clinically proven to reduce colic

Care

- Moisturising nipple cream
- Breast shell set
- Unique four layered construction for maximum dryness

Highlights

Electronic breast pump



Only the Philips Avent electric breast pump has a unique memory feature that learns and continues your personal pumping rhythm for maximum comfort. The breast pump also has unique active massage cushions, which are designed to help stimulate natural let-down.

Bottles clinically proven



Unique anti-colic system consists of the Philips Avent teat and adaptor ring allowing air to flow into the bottle instead of your baby's tummy for better digestion.*

Milk storage cups



Ideal for storing breast milk: Thanks to the adaptor that is compatible with Philips Avent breast pumps and teats, you can express your milk directly into the cup and feed your baby from the same cup. Cups are pre-sterilised/ready to use, easy to label, leak proof, can be used in a fridge/freezer and are dishwasher safe

Moisturising nipple cream



Hypoallergenic for mother and baby. Moisturises and softens dry or sensitive nipples. Easy to apply and no need to remove before breastfeeding.

Breast shell set

The ultra-soft breast shells are worn inside your bra to protect your nipples from chafing and to collect excess breast milk

Disposable breast pads



Unique four layered construction for maximum dryness: 1. Ultra-soft top layer with nipple indent — keeps breast dry at all times 2. Ultra-absorbent layer, cushioned for extra comfort 3. Ultra-absorbent core — draws in moisture to guard against embarrassment 4. Ultra-breathable outer layer — helps avoid sore nipples. Anti-slip with adhesive tape — keeps pad in place

Sealing discs



Convert all Philips Avent bottles into milk storage containers.

Specifications

Country of origin

EU

What is included

Electronic Breast Pump: 1

Classic Bottle 9 oz: 2 pcs

Classic Bottle 4 oz: 1 pcs

Extra-soft slow flow teat: 2 pcs

Newborn teat travel pack: 1 pcs

Storage Cup (180 ml/6 oz): 10 pcs

Milk storage cup adapter: 2 pcs

Milk storage cup lid: 10 pcs

Sealing disc for milk storage: 1 pcs

Moisturising Nipple Cream: 1 pcs

Breast Shell: 6 pcs

Disposable breast pad: 40 pcs

Guide to breastfeeding DVD: 1 pcs

Development stages

Stages: 0–6 months

