

# PHILIPS

## AVENT

Baby feeding  
essentials

SCD285/00



## Baby feeding made easy

Express, store, feed, clean, sterilise and keep your bottles cold or warm when travelling, thanks to our Baby Feeding Essentials Set.

### Clean & Sterilise

- Curved brush head and moulded tip
- Microwave steam steriliser

### Express

- Manual breast pump with cup

### Feed

- Bottles clinically proven to reduce colic\*

### Store

- Milk storage cups

### Travel

- Neoprene ThermoBag (Black)
- Milk powder dispenser

# Highlights

## Microwave steam steriliser



Compact and lightweight (fits most microwaves); convenient for travel; large capacity (holds 6 Philips AVENT bottles); ultra fast; easy to use (just add water, load and place in the microwave); safe, easy handling (side grips close the lid securely)

## Manual breast pump with cup



Our manual breast pump has a unique active massage cushion intended to help stimulate natural let-down. It has been designed to increase your comfort when you are pumping. Express in a bottle or directly into a storage cup. Thanks to the adaptor compatible with Philips AVENT teats, you can also feed your baby directly from the cup, no need to transfer into a bottle.

## Bottle and teat brush



Specially curved brush head and moulded handle-tip reach the corners of all types of bottles, teats and feeding equipment for thorough cleaning.

## Milk storage cups



Ideal for storing breast milk: Thanks to the adaptor that is compatible with Philips AVENT breast pumps and teats, you can express your milk directly into the cup and feed your baby from the same cup. Cups are pre-sterilised/ready to use, easy to label, leak proof, can be used in a fridge/freezer and are dishwasher safe

## Bottles clinically proven



Unique anti-colic system consists of the Philips AVENT teat and adaptor ring allowing air to flow into the bottle instead of your baby's tummy for better digestion.

## Milk powder dispenser



This handy little unit carries 3 pre-measured portions of milk powder in separate compartments. When you're ready to feed, just pour the powder into the bottle of pre-boiled cooled water. Remove inner sections to use as a bowl or container.

## Neoprene ThermoBag (Black)



Holds 2 Philips AVENT Bottles, 2 toddler cups or 4 storage cups and ensures bottles stay warm or cool (storage of cold expressed breast milk/formula or pre-boiled hot water for up to 4 hours).

# Specifications

## What is included

Microwave steam steriliser: 1 pcs  
Manual breast pump + cup and lid: 1 pcs  
Classic Bottle 4 oz: 2 pcs  
Classic Bottle 9 oz: 4 pcs  
Extra soft Newborn Flow Teat: 2 pcs  
Extra soft Slow Flow Teat: 4 pcs  
Soft, Fast Flow teat: 2 pcs

Newborn teat travel pack: 1 pcs  
Storage Cup (180 ml/6 oz): 1 pcs  
Milk storage cup adapter: 1 pcs  
Lid: 1 pcs  
ThermaBag: 1 pcs  
Milk powder dispenser: 1 pcs  
Bottle and teat brush: 1 pcs

## Country of origin

EU

## Development stages

Stages: 0–6 months

