

# PHILIPS

## AVENT

Newborn Starter Set

Natural

SCD290/00



## The natural way to start bottle feeding

### Bottle, teat and soother set

A handy collection including 4 Natural bottles (2 x 4 oz and 2 x 9 oz), a bottle and teat brush, and a white Advanced orthodontic soother 0-6 months (designed to help healthy oral development). The new bottle makes bottle feeding more natural

#### Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast-shaped teat

#### A more comfortable and contented feed for baby

- Unique Petals for a soft, flexible teat without collapse

#### Vents air away from your baby's tummy

- Advanced anti-colic system with innovative twin valve

#### Easy to hold and grip in any direction

- Ergonomic shape for maximum comfort

#### Other benefits

- Advanced orthodontic soother
- Curved brush head for easy cleaning
- Compatible with the Philips Avent Natural range
- Simple to use and clean, quick and easy assembly
- These bottles are BPA free\*



[asimpleswitch.com](http://asimpleswitch.com)

# Highlights

## Natural latch on



The wide, breast-shaped teat promotes natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

## Ergonomic shape




Due to the unique shape, the bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

## Curved brush head



Specially curved brush head and moulded handle-tip reach the corners of all wide neck bottles, teats and feeding products for thorough cleaning.

## Unique comfort petals



Petals inside the teat increase softness and flexibility without teat collapse. Your baby will enjoy a more comfortable and contented feed.

## Advanced orthodontic soother



Developed with leading orthodontist Dr. Hagemann, the advanced orthodontic soother is designed to help healthy oral development. It features wings specially shaped to minimise pressure on gums and developing teeth.

## Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

## Advanced anti-colic valve



Innovative twin valve designed to reduce colic and discomfort by venting air into the bottle and not baby's tummy.

## Compatibility across the range



The new Philips Avent Natural bottle is compatible with the Philips Avent range, excluding Classic bottles and cup handles. We advise using the Natural bottles with Natural feeding teats only.



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Country of origin

England

## Material

**Bottle:** Polypropylene, BPA free\*

**Teat:** Silicone, BPA free\*

## What is included

**Baby bottle:** 4 pcs

**Bottle and teat brush:** 1 pcs

**Advanced orthodontic soother:** 1 pcs

## Bottle

**Material:** BPA free\*

## Design

**Bottle design:** Ergonomic shape, Wide neck

## Ease of use

**Bottle use:** Easy to hold, Easy to assemble,  
Easy to clean

## Development stages

**Stages:** 0–6 months

## Functions

**Latch on:** Easily combine breast and bottle,  
Natural latch on

**Teat:** Extra-soft and flexible teat, Unique  
comfort petals

**Anti-colic valve:** Advanced anti-colic system

## Features

**Can be sterilised**

