

# PHILIPS

## AVENT

Newborn Natural  
starter set

Natural

Newborn Starter Set



SCD290/11



## The most natural way to bottle feed

### Natural latch on

Our new bottle with the skin soft nipple material and flexible spiral design, more closely resembles the breast. The comfort petals and natural nipple shape allows natural latch on and makes it easy to combine breast and bottle feeding.

#### Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple
- Skin-soft teat designed for a natural feed

#### A comfortable and contented feed for your baby

- Flexible spiral design combined with comfort petals

#### Vents air away from your baby's tummy

- Unique anti-colic valve technology

#### Other benefits

- Designed for your baby's changing needs
- Ergonomic shape for maximum comfort
- Compatible with the Philips Avent Natural range
- This bottle is BPA free\*
- Simple to use and clean, quick and easy assembly
- Orthodontic pacifier
- Bottle brush with curved brush head for easy cleaning

# Highlights

## Natural latch on



The wide breast shaped nipple promotes a natural latch on similar to the breast. Designed to give a more comfortable and contented feed for your baby.

## Skin-soft teat



The skin-soft nipple material closely resembles the feeling of the breast, making it easy for your baby to combine breast and bottle feeding.

## Flexible spiral teat design



The spiral design combined with our comfort petals inside the nipple increases softness and the flexibility allowing natural tongue movement without nipple collapse. Designed to give your baby a more comfortable and content feeding.

## Anti-colic valve technology



Reduces fussing and discomfort by venting air away from baby's tummy.

## Baby's changing needs

The Philips Avent Natural bottle is available in 4 sizes and 7 different nipples for every development stage of your baby. With different nipple softness, flexibility and design and with increasing flow rates and bottle sizes that keep up with your baby's growth and development.

## Ergonomic shape



Due to the unique shape, the bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

## Compatibility across the range



The new Philips Avent Natural bottle is compatible with the Philips Avent range, excluding Classic bottles and cup handles. We advise to use the Natural bottles with Natural feeding nipples only.

## This bottle is BPA free\*



The Philips Avent Natural bottle is made of BPA free\* material (polypropylene).

## Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

## Orthodontic pacifier



The orthodontic, collapsible nipple of the translucent pacifier respects the natural development of baby's palate, teeth and gums. All our pacifiers are made from silicone and are taste and odour-free.

# Specifications

## Material

**Bottle:** Polypropylene, BPA free\*

**Nipple:** Silicone, BPA free\*

## What is included

**4oz Baby bottles:** 2 pcs

**9oz Baby bottle:** 2 pcs

**Soother:** 1 pcs

**Brush:** 1 pcs

## Bottle

**Material:** BPA free\*

## Design

**Bottle design:** Ergonomic shape, Wide neck

## Ease of use

**Bottle use:** Easy to hold, Easy to assemble,  
Easy to clean

## Development stages

**Stage:** 0-6 months

## Functions

**Latch on:** Natural latch on, Easy combine  
breast and bottle

**Nipple:** Unique comfort petals, Extra soft and  
flexible teat

**Anti-colic valve:** Advanced anti-colic system

