

PHILIPS

AVENT

Newborn Glass
Starter Set

Natural

SCD291/01



The most natural way to bottle feed

Easy to combine with breastfeeding

The Philips Avent Natural glass newborn starter set has all you need to start bottle feeding in a natural way. It includes 5 Natural glass bottles in multiple sizes (4oz and 8oz), a glass bottle sleeve and 2 soothie pacifiers.

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple
- Unique Petals for a soft, flexible nipple without collapse

Vents air away from your baby's tummy

- Advanced anti-colic system with innovative twin valve

Other benefits

- Pharmaceutical grade glass
- Compatible with the Philips Avent Natural range
- Easily read the scale of the glass bottle
- Heat resistance
- This bottle is BPA free*

Other essential accessories included in set

- Bottle sleeve protects 4oz/120ml glass bottles
- One piece, silicone Soothie helps calm babies

Highlights

Natural latch on



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Unique comfort petals



Petals inside the nipple increase softness and flexibility without nipple collapse. Your baby will enjoy a more comfortable and contented feed.

Advanced anti-colic valve



Innovative twin valve designed to reduce colic and discomfort by venting air into the bottle and not baby's tummy.

Pharmaceutical grade glass

Premium quality borosilicate glass to ensure the highest quality for ultimate purity.

Compatibility across the range



The new Philips Avent Natural bottle is compatible with the Philips Avent range, excluding Classic bottles and cup handles. We advise to use the Natural bottles with Natural feeding nipples only.

Easily read the scale

The opening in the Philips Avent glass bottle sleeve, means you can easily read the scale during your baby's feed.

Heat resistance



The Natural glass bottle is heat and thermal shock resistant. It can therefore be safely stored in the fridge, warmed, and is also suitable for sterilization.

This bottle is BPA free*



The Philips Avent Natural bottle is made of BPA free* material (glass, polypropylene and silicone).

Protective bottle sleeve



The insulating bottle sleeve, allows you to keep the feed at the right temperature and to protect your glass bottle at the same time.

Silicone soothie



The silicone soothie helps calm your baby. It is used by medical professionals to calm newborns.

Specifications

Design

Bottle design: Ergonomic shape, Wide neck

Development stages

Stage: 0 - 6 months

Ease of use

Bottle: Dishwasher & microwave safe, Easy to assemble, Easy to clean, Easy to hold

Easy to use

Sleeve use: Easy scale reading, Fits glass bottles, Shock resistant, Warm and cool insulating

Functions

Anti-colic valve: Advanced anti-colic system

Latch on: Natural latch on, Easy combine breast and bottle

Nipple: Unique comfort petals, Extra soft and flexible nipple

Material

Bottle: BPA free*, Glass

BPA free*: Bottle made of Borosilicate type 1 glass; Nipple made of silicone; Screw ring made of PP (Polypropylene)

Nipple: BPA free*, Silicone

Sleeve: Oeko-Tex® Class 1 certified, Soft polyester

What is included

Baby Bottle: 5 pcs

Glass bottle sleeve: 1 pcs

Soothie pacifier: 2 pcs

