

PHILIPS

All in One

Natural

SCD293/00



Sterilise, Express and Feed

Breastfeeding made easy

The Natural All-in-one starter set is the ideal set for mums who want to start expressing and feeding in a clean and simple way. The set contains all the products you need to sterilise, express your milk and feed your baby.

Express

- More comfortable pumping position due to unique design
- Soft massage cushion with warm feel

Effective sterilisation

- Natural steam sterilisation kills 99.9% of harmful germs

3-in-1 steriliser

- Adjustable size

Fits all bottle types

- Suitable for wide neck and narrow neck bottles

Feed

- Natural latch on due to the wide breast-shaped teat
- Unique Petals for a soft, flexible teat without collapse

Highlights

Natural latch on



The wide, breast-shaped teat promotes natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means that you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

Soft massage cushion



Our massage cushion has a new soft, velvety texture that gives the skin a warm feeling for comfortable, gentle stimulation of your milk

flow. The cushion is designed to gently mimic your baby's suckling to help stimulate let-down.

Fits all bottle types



The steriliser has been designed to fit both wide neck and narrow neck bottles as well as different breast pumps and accessories. This gives you the freedom to choose the bottles and accessories that best suit you and your baby. At full capacity the steriliser holds six 11 oz/330 ml Philips Avent bottles.

3-in-1 steriliser



Due to its adjustable size, this 3-in-1 Philips Avent steriliser perfectly fits the items you would like to sterilise and always takes up the smallest amount of kitchen space. The 3-in-1 functionality gives you: 1) Small size configuration for sterilising soothers 2) Medium size configuration for sterilising breast pumps, toddler plates, knives and forks 3) Large size configuration for sterilising 6 baby

Kills 99.9% of harmful germs



Sterilising protects your baby from particularly harmful bacteria until his immune system is strong enough. Philips Avent Sterilisers use the hospital method of steam sterilisation, which is quick, easy and effective, with no chemicals involved. Contents stay sterile for 24 hours if the lid is unopened.

Unique comfort petals



Petals inside the teat increase softness and flexibility without teat collapse. Your baby will enjoy a more comfortable and contented feed.

Specifications

Country of origin

England
Turkey

Material

Bottle: Polypropylene, BPA free*
Breast pump: BPA free*, Polypropylene

What is included

3-in-1 steriliser: 1 pcs
Bottle and teat brush: 1 pcs
Extra soft Newborn Flow teat: 2 pcs
Sealing disc for milk storage: 1 pcs

Disposable day breast pads: 2 pcs
Disposable night breast pads: 2 pcs
Natural bottle 125 ml/4 oz: 2 pcs
Natural bottle 9 oz: 2 pcs
Silicone Soother: 3 pcs
Snap on hygienic cap: 3 pcs
Tongs: 1 pcs
Breast pump body with handle: 1 pcs

Development stages

Stages: 0–6 months

Functions

No leaning forward: Sit in a comfortable position
Soft massage cushion: Gentle stimulation

Power

Voltage: 220 V V

Compatibility

Philips Avent range-compatible

