

PHILIPS

AVENT

Newborn Natural
starter set

Newborn Starter Set



SCD301/01



Easy to combine with breastfeeding Natural latch on

Our Natural bottle with an ultra-soft teat more closely resembles the breast. The wide breast-shaped teat with flexible spiral design and comfort petals allows natural latch-on and makes it easy to combine breast and bottle feeding.

Easy to combine breast and bottle feeding

- Natural latch-on due to the wide breast shaped teat
- Ultra soft teat designed to mimic the feel of the breast

A comfortable and contented feed for your baby

- Flexible spiral design combined with comfort petals

Vents air away from your baby's tummy

- Anti-colic valve designed to reduce colic and discomfort

Other benefits

- Teats with different flow rates available
- Ergonomic shape for maximum comfort
- Compatible with the Philips Avent Natural range
- This bottle is BPA free*
- Simple to use and clean, quick and easy assembly
- Orthodontic soother
- Bottle brush with curved brush head for easy cleaning

Highlights

Natural latch on



The wide breast shaped teat promotes a natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Anti-colic valve



Anti-colic valve designed to keep air away from your baby's tummy, to help reduce colic and discomfort.

This bottle is BPA free*



The Philips Avent Natural bottle is made of BPA-free* material (polypropylene).

Ultra soft teat



The teat has an ultra soft texture, designed to mimic the feel of the breast.

Choosing the right teat



The Philips Avent Natural range offers different teat softness and increasing flow rates for each of your baby's development stages.

Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

Flexible spiral design



Flexible spiral design, combined with our unique comfort petals to create a flexible teat, allowing for a more natural feed without teat collapse.

Ergonomic shape



Due to its unique shape, the bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

Orthodontic soother



The orthodontic, collapsible teat of the soother respects the natural development of baby's palate, teeth and gums. All our soothers are made from silicone and are taste and odour-free.

Compatibility across the range



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

Specifications

Country of origin

United Kingdom

Material

Bottle: Polypropylene, BPA-free*

Teat: Silicone, BPA-free*

What is included

4 oz Baby bottles: 2 pcs

9 oz Baby bottle: 2 pcs

Brush: 1 pcs

Soother: 1

Design

Bottle design: Ergonomic shape, Wide neck

Ease of use

Bottle use: Easy to hold, Easy to assemble,

Easy to clean

Development stages

Stages: 0-12 months

Functions

Latch on: Natural latch on, Easily combine breast and bottle

Teat: Unique comfort petals, Ultra soft and flexible teat

Anti-colic valve



* 0% BPA, following EU regulation 10/2011

* What is colic and how does it affect babies? Colic is caused in part by swallowing air while feeding, which creates discomfort in a baby's digestive system. Symptoms include crying and fussing.