



Newborn Natural starter set

- Newborn Starter Set

SCD301/01

Easy to combine with breastfeeding

Natural latch on

Our Natural bottle with an ultra-soft teat more closely resembles the breast. The wide breast-shaped teat with flexible spiral design and comfort petals allows natural latch-on and makes it easy to combine breast and bottle feeding.

Benefits

Easy to combine breast and bottle feeding

- Natural latch-on due to the wide breast shaped teat
- Ultra soft teat designed to mimic the feel of the breast

A comfortable and contented feed for your baby

- Flexible spiral design combined with comfort petals

Vents air away from your baby's tummy

- Anti-colic valve designed to reduce colic and discomfort

Other benefits

- Teats with different flow rates available
- Ergonomic shape for maximum comfort
- Compatible with the Philips Avent Natural range
- This bottle is BPA free*
- Simple to use and clean, quick and easy assembly
- Orthodontic soother
- Bottle brush with curved brush head for easy cleaning

PHILIPS

Features

Natural latch on



The wide breast shaped teat promotes a natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Ultra soft teat



The teat has an ultra soft texture, designed to mimic the feel of the breast.

Flexible spiral design



Flexible spiral design, combined with our unique comfort petals to create a flexible teat, allowing for a more natural feed without teat collapse.

Anti-colic valve



Anti-colic valve designed to keep air away from your baby's tummy, to help reduce colic and discomfort.

Choosing the right teat



The Philips Avent Natural range offers different teat softness and increasing flow rates for each of your baby's development stages.

Ergonomic shape



Due to its unique shape, the bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

Compatibility across the range



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

This bottle is BPA free*



The Philips Avent Natural bottle is made of BPA-free* material (polypropylene).

Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

Orthodontic soother



The orthodontic, collapsible teat of the soother respects the natural development of baby's palate, teeth and gums. All our soothers are made from silicone and are taste and odour-free.

Curved brush head



Specially curved brush head and moulded handle-tip reach the corners of all wide neck bottles, teats and feeding products for thorough cleaning.

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Country of origin

United Kingdom Yes

Material

Bottle Polypropylene
BPA-free*
Teat Silicone
BPA-free*

What is included

4 oz Baby bottles 2 pcs
9 oz Baby bottle 2 pcs
Brush 1 pcs
Soother 1

Design

Bottle design Ergonomic shape
Wide neck

Ease of use

Bottle use Easy to hold
Easy to assemble
Easy to clean

Development stages

Stages 0-12 months

Functions

Latch on Natural latch on
Easily combine breast and bottle

Teat Unique comfort petals, Ultra soft and flexible teat
Anti-colic valve Yes

Packaging dimensions

Height 18.30 cm
Width 28.20 cm
Depth 7.20 cm
Net weight 0.39 kg
Gross weight 0.46 kg
EAN 08710103868880
Number of products included 1
Country of origin GB
Harmonised Systems Code 392410

Outer Carton

Length 36.40 cm
Width 29.10 cm
Height 19.50 cm
Gross weight 2.53 kg
EAN 18710103868887
Number of consumer packages 5

* 0% BPA, following EU regulation 10/2011

* What is colic and how does it affect babies? Colic is caused in part by swallowing air while feeding, which creates discomfort in a baby's digestive system. Symptoms include crying and fussing.



data subject to change
2019, July 10

Version: 3.0.1
EAN: 08710103868880

© 2019 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com