

PHILIPS

AVENT

Newborn Starter Set

Pink

SCD363/02



Helps settle your baby, especially at night*

Less colic*

Philips Avent Classic+ Newborn Starter Set is the ideal set to get started with bottle feeding. It includes 4 Classic+ feeding bottles (2x 4Oz, 2x 9Oz), brush, milk powder dispenser, 2 pacifiers and a special trainer solution.

Benefits and features:

- Curved brush head and molded tip
- Bottles clinically proven to reduce colic
- Holds enough powdered formula for three 260 ml/ 9 oz feeds
- Orthodontic, symmetrical and collapsible nipple

0% BPA (Bisphenol A Free)

- This bottle is BPA free

Vents air away from your baby's tummy

- Clinically proven anti-colic system

More comfortable feed at baby's chosen pace

- Easy latch on due to the unique valve on the nipple

Highlights

Bottle and nipple brush




Specially curved brush head and molded handle-tip reach the corners of all types of bottles, nipples and feeding equipment for thorough cleaning.

Easy latch on



Unique valve on the nipple flexes to your baby's feeding rhythm. Milk will only flow at the pace chosen by your baby to minimize overeating and spit-up, burping and gas

Newborn Silicone Pacifier



Flat, drop-shaped symmetrical nipples respect the natural development of your baby's palate, teeth and gums, even if the pacifier ends up upside-down in the mouth.

Bottles clinically proven



Unique anti-colic system consists of the Philips Avent nipple and adaptor ring allowing air to flow into the bottle instead of your baby's tummy for better digestion.*

Holds enough formula for three 260 ml/ 9 oz feeds



The Philips Avent milk powder dispenser holds 3 pre-measured doses of milk powder - ideal for travel

Unique anti-colic system



As your baby feeds, the unique valve on the nipple flexes to allow air into the bottle instead of your baby's tummy.*

This bottle is BPA free



The new Philips Avent Natural bottle is made of BPA free material (polypropylene)

Specifications

Country of origin

England

What is included

Bottle and nipple brush: 1 pcs

Milk powder dispenser: 1 pcs

Baby bottle: 4 pcs

Trainer cup: 1 pcs

Spout: 1 pcs

Pacifier: 2 pcs

Bottle

Material: BPA free

Design

Bottle design: Wide neck, Ergonomic shape

Ease of use

Bottle use: Easy to clean, Easy to hold

Development stages

Stages: 0-6 months

Functions

Anti-colic valve: Integrated anti-colic system

Latch on: Easy latch on, Promotes nipple acceptance

Nipple: Flexes to feeding rhythm

