

PHILIPS

Baby gift set

Infant Starter Set



SCD364/01



Clinically proven to reduce colic and discomfort* Easy to clean, healthier feeding

Our trusted Classic range has been improved to provide an even more enjoyable feeding experience. The clinically proven, anti-colic system is now integrated into the nipple, making the bottle extra easy to clean and assemble.

Clinically proven to reduce colic and discomfort*

- Anti-colic system proven to reduce colic*
- Fussing is reduced, especially at night
- Baby controls milk flow for less spit-up, burping and gas

Benefits and features:

- This bottle is BPA free
- Curved brush head and molded tip
- Holds enough powdered formula for three 260 ml/ 9 oz feeds
- Orthodontic pacifier

Highlights

Less colic



Unlike other bottles, the clinically proven anti-colic system is now integrated into the nipple, making it easier to assemble the bottle correctly. As your baby feeds, the unique valve on the nipple opens to allow air into the bottle instead of your baby's tummy.

Less fussing*



Sleep and nutrition are vital to your baby's health and happiness. A randomised clinical trial was carried out to see whether infant baby bottle design affects "infant behavior". The Philips Avent Classic baby bottle was shown to significantly reduce fussing by approximately 28 minutes a day as compared to the comparator bottle (46 min vs 74 min, p=0.05). This was especially true during the night-time.*

Less spit-up, burping and gas




The unique valve on the nipple flexes to your baby's feeding rhythm. Milk will only flow at the pace chosen by your baby to minimise overeating and spit-up, burping and gas.

BPA free



Philips Avent anti-colic bottle is made of BPA free material (PP).

Bottle and nipple brush



Specially curved brush head and molded handle-tip reach the corners of all types of bottles, nipples and feeding equipment for thorough cleaning.

Holds enough formula for three 260 ml/ 9 oz feeds



The Philips Avent milk powder dispenser holds 3 pre-measured doses of milk powder - ideal for travel

Orthodontic pacifier



The orthodontic, collapsible nipple of the translucent pacifier respects the natural development of baby's palate, teeth and gums. All our pacifiers are made from silicone and are taste and odour-free.

Specifications

Bottle

Material: BPA free*

Design

Bottle design: Ergonomic shape, Wide neck

Development stages

Stages: 0-12 months

Ease of use

Ease of use: Dishwasher & microwave safe,
Easy to assemble, Easy to clean, Easy to hold

Functions

Clinically proven: clinically proven anti-colic

Nipple: Easy latch on, Flexes to feeding
rhythm

Material

Bottle: BPA free*, Polypropylene

Nipple: BPA free*, Silicone

What is included

Bottle and nipple brush: 1 pcs

Milk powder dispenser: 1 pcs

Translucent pacifier: 1 pcs

Trainer cup: 1 pcs

Baby Bottle: 4 pcs

