

PHILIPS

AVENT

Newborn Starter Set

Classic+ (new)

SCD373



Clinically proven to reduce colic and discomfort*

Easy to clean for perfect hygiene

The handy collection SCD373/00 by Philips Avent including 5 Classic+ feeding bottles (2x125 ml, 2x260, 1x330 ml), a bottle and nipple brush, and a white translucent pacifier 0–6 months.

Benefits and features:

- Curved brush head and molded tip

Other benefits

- BPA free
- Compatible with the Philips Avent range
- Different nipple flow rates available

Clinically proven to reduce colic and discomfort*

- Anti-colic system proven to reduce colic*
- Fussing is reduced, especially at night
- Baby controls milk flow for less spit-up, burping and gas

Easy to use, clean and assemble

- Ergonomic shape for easy holding
- Few parts for a quick and simple assembly
- Wide neck and only few parts for easy cleaning
- Leak-free for an enjoyable feeding


Highlights

Less colic



Unlike other bottles, the clinically proven anti-colic system is now integrated into the nipple, making it easier to assemble the bottle correctly. As your baby feeds, the unique valve on the nipple opens to allow air into the bottle instead of your baby's tummy.

Less fussing*



Sleep and nutrition are vital to your baby's health and happiness. A randomised clinical trial was carried out to see whether infant baby bottle design affects "infant behavior". The Philips Avent Classic baby bottle was shown to significantly reduce fussing by approximately 28 minutes a day as compared to the comparator bottle (46 min vs 74 min, $p=0.05$). This was especially true during the night-time.*

Less spit-up, burping and gas



The unique valve on the nipple flexes to your baby's feeding rhythm. Milk will only flow at the pace chosen by your baby to minimise overeating and spit-up, burping and gas.

Leak-free



The new Classic+ bottle is designed to prevent leakage whilst feeding, for a truly enjoyable feeding experience.

Easy to clean



With only 4 parts, a wide bottle neck, and rounded corners, our bottle is easy to clean quickly and thoroughly. Have peace of mind that your bottle is truly clean in no time.

Easy to assemble in seconds



The new Classic+ bottle has only 4 parts for quick and simple assembly.

Easy to hold



Due to the unique shape, the feeding bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

BPA free

Following EU directive (2011/8/EU), the Philips Avent Classic+ feeding bottle is made of BPA free material (polypropylene).

Compatibility across range



The Philips Avent Classic+ bottle is compatible with most of the Philips Avent range. We recommend to use the Classic+ bottle with Classic+ teats only.

Fits your baby's growing needs



The Philips Avent Classic+ range offers different nipple flow rates to keep up with your baby's growth. Remember that age indications are approximate as babies develop at different rates. All nipples are available in twin packs: Newborn, Slow, Medium, Fast and Variable Flow.

Bottle and nipple brush



Specially curved brush head and molded handle-tip reach the corners of all types of bottles, nipples and feeding equipment for thorough cleaning.

Specifications

Country of origin

England

Material

Bottle: Polypropylene, BPA free*

Nipple: Silicone, BPA free*

What is included

Bottle and nipple brush: 1 pcs

Translucent pacifier: 1 pcs

Baby Bottle: 5 pcs

Bottle

Material: BPA free*

Design

Bottle design: Wide neck, Ergonomic shape

Ease of use

Bottle use: Easy to clean, Easy to hold

Ease of use: Easy to hold, Dishwasher & microwave safe, Easy to assemble, Easy to clean

Development stages

Stages: 0-12 months

Functions

Clinically proven: clinically proven anti-colic

Nipple: Flexes to feeding rhythm, Easy latch on

Convenience

Capacity: 4, 9, 11 oz

Features

Can be sterilized

