


PHILIPS

AVENT

Bottle feeding
essentials

Classic+ (new)

SCD383/01



Bottle feeding made easy

Easy to clean for perfect hygiene

Clean and sterilise your bottles and teats, feed your baby and keep your bottles cold or warm when travelling, thanks to the Philips Avent Bottle Feeding Essentials Set.

Feed

- Bottles clinically proven to reduce colic

Other benefits

- BPA free

Clinically proven to reduce colic and discomfort*

- Anti-colic system proven to reduce colic*
- Fussing is reduced, especially at night
- Baby controls milk flow for less spit-up, burping and wind

Easy to use, clean and assemble

- Ergonomic shape for easy holding
- Few parts for quick and simple assembly
- Wide neck and only a few parts for easy cleaning
- Leak-free for enjoyable feeding

Highlights

Bottles clinically proven



Unique anti-colic system consists of the Philips Avent teat and adaptor ring allowing air to flow into the bottle instead of your baby's tummy for better digestion.*

Easy to hold



Due to the unique shape, the feeding bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

Less colic



In contrast to other bottles, the clinically proven anti-colic system is now integrated into the teat, making it easier to assemble the bottle correctly. As your baby feeds, the unique valve on the teat opens to allow air into the bottle instead of your baby's tummy.

BPA free

In accordance with EU directive 2011/8/EU, the Philips Avent Classic+ feeding bottle is made of BPA-free material (polypropylene).

Easy to assemble in seconds




The new Classic+ bottle has only 4 parts for quick and simple assembly.

Easy to clean



With only 4 parts, a wide bottle neck and rounded corners, our bottle is easy to clean quickly and thoroughly. Have peace of mind that your bottle is truly clean in no time.

Less fussing*



Sleep and nutrition are vital to your baby's health and happiness. A randomised clinical trial was carried out to see whether infant baby

bottle design affects "infant behaviour". The Philips Avent Classic baby bottle was shown to significantly reduce fussing by approximately 28 minutes a day as compared to the comparator bottle (46 min vs. 74 min, p=0.05) This was especially true at night.*

Leak-free



The new Classic+ bottle is designed to prevent leakage whilst feeding, for a truly enjoyable feeding experience.

Less spit-up, burping and wind



The unique valve on the teat flexes to your baby's feeding rhythm. Milk will only flow at the pace chosen by your baby to minimise overeating and spit-up, burping and wind.

Specifications

Material

Bottle: Polypropylene, BPA-free*

Teat: Silicone, BPA-free*

What is included

Extra-soft newborn flow teat: 2 pcs

Bottle and teat brush: 1 pcs

Extra-soft slow flow teat: 4 pcs

Milk powder dispenser: 1 pcs

Soft, Fast Flow teat: 2 pcs

Microwave steam steriliser: 1 pcs

ThermaBag: 1 pcs

Baby Bottle: 6 pcs

Bottle

Material: BPA-free*

Design

Bottle design: Wide neck, Ergonomic shape

Ease of use

Bottle use: Easy to clean, Easy to hold

Ease of use: Dishwasher- and microwave-safe, Easy to assemble, Easy to clean, Easy to hold

Development stages

Stage: 0–6 months, 0–12 months

Functions

Clinically proven: clinically proven anti-colic

Teat: Flexes to feeding rhythm, Easy latch on

Features

Can be sterilised

