

PHILIPS

AVENT

Anti-colic Bottle with
AirFree vent Gift Set

Newborn Starter Set

Bottles with AirFree vent

Pink

SCD393/04



Designed to reduce colic, gas and reflux*

Our unique AirFree vent is designed to help your baby swallow less air as they drink by keeping the nipple full of milk during feeding. Reducing the amount of air your baby ingests helps w/common feeding issues such as colic, reflux & gas.

AirFree vent designed for less air ingestion

- Nipple stays full of milk, not air
- Helps baby swallow less air with AirFree vent

Additional benefits

- Easy to clean - wide neck & few parts
- AirFree vent is one piece, simple to use
- Use your bottle with or without the AirFree vent
- Different nipple flow rates available
- Mix and match with other Philips Avent products
- Dispenser holds enough formula for three 9oz/260ml feeds
- Orthodontic, symmetrical and collapsible pacifier nipple

Clinically proven to reduce colic and fussing*


- Valve in the nipple prevents vacuum build-up*
- Nipple shape designed for a secure latch



asimpleswitch.com


Highlights

Nipple full of milk, not air



The nipple stays full even when the bottle is horizontal, so your baby can drink in an upright position. This can help reduce reflux, aid digestion, and make feeding time more comfortable for you and your little one.

Reduces colic and fussing*



Clinical studies have shown that the Philips Avent bottle reduces colic and fussing. Fussing is significantly reduced at night as babies fed with Philips Avent anti-colic bottles experienced 60% less fussing than babies fed with a leading competitor's vented bottle.

Designed for a secure latch



The shape of the nipple allows for a secure latch, and the ribbed texture helps prevent nipple collapse, for uninterrupted, comfortable feeding.

Unique venting system




Our unique AirFree vent is designed to help your baby swallow less air as they drink, by drawing air away from the nipple. Reducing the amount of air your baby ingests can help with common feeding issues, such as colic, reflux and gas.

Simple to use



While feeding, make sure the AirFree vent stays on top. This keeps the nipple full of milk, even when the bottle is held in a horizontal position, for easier upright feeding.

Choosing the right nipple



Different flow rates to keep up with your baby's growth. The Philips Avent Anti-colic bottles can be used with or without the AirFree vent – but they should only be used with Philips Avent Anti-colic feeding nipples.

Easy to clean



With its wide neck and fewer parts, our bottle is easier to assemble and to clean quickly and thoroughly.

Freedom of choice



You can use your Philips Avent Anti-colic bottle with or without the AirFree vent. The AirFree vent is compatible with all of our Philips Avent Anti-colic bottles.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Country of origin

Assembled in the USA

England: Bottle(s), Pacifier(s)

Malaysia: Bottle brush, Formula dispenser

Indonesia: Spout

Material

Bottle: Polypropylene, BPA free

Nipple: Silicone, BPA free

What is included

Baby bottle: 4 pcs

Bottle and nipple brush: 1 pcs

Milk powder dispenser: 1 pcs

Trainer cup: 1 pcs

Spout: 1 pcs

Pacifier: 2 pcs

AirFree vent: 4 pcs

Design

Bottle design: Ergonomic shape, Wide neck

Ease of use

Bottle use: Dishwasher & microwave safe, Easy to assemble, Easy to clean, Easy to hold

Development stages

Stages: 0-12 months

Functions

Nipple: Easy latch on, Ribbed texture prevents nipple collapse, Proven Anti-colic system

Ease of use: Easy to clean and assemble, Few parts for easy assembly

Compatibility: Full Classic+ compatibility



- * Our unique AirFree vent is designed to help your baby swallow less air as they drink by keeping the nipple full of milk during feeding. Reducing the amount of air your baby ingests helps w/common feeding issues such as colic, reflux & gas.
- * At 2 weeks of age, babies fed with a Philips Avent bottle showed a trend to less colic compared to a conventional bottle and a significant reduction in fussing at night compared to babies fed with another leading bottle.
- * What colic is, and how it affects babies? Colic is caused in part by swallowing air while feeding, which creates discomfort in a baby's digestive system. Symptoms include crying, fussing, gassiness and spit-up.