

# PHILIPS

## AVENT

Anti-colic bottle gift set

Infant Starter Set

SCD394/01



## Clinically proven to reduce colic and discomfort\*

Designed to reduce gas

Our anti-colic range's Airflex venting system and textured nipple is designed to minimize feeding interruptions and discomfort. With its integrated anti-colic valve, air is vented into the bottle and away from the baby's tummy.

### **Clinically proven to reduce colic and discomfort\***

- Anti-colic system proven to reduce colic\*
- 60% less fussing at night\*
- Ribbed texture prevents collapse for uninterrupted feeding
- Nipple shape designed for secure latch

### **Easy to use, clean and assemble**

- Leak-free design
- Easy to assemble with only fewer parts
- Easy to clean bottles especially with curved brush head

### **Other benefits**

- Compatible range from breastfeeding to cup
- Orthodontic, symmetrical and collapsible pacifier nipple
- This bottle is BPA free
- Full compatibility with Classic+ bottle range
- Dispenser holds enough formula for three 260ml/9oz feeds


# Highlights

## Less colic



Our clinically proven Airflex venting system is designed to reduce gas and colic. As your baby feeds, the valve integrated into the nipple flexes to allow air into the bottle to prevent vacuum build up and vents it towards the back of the bottle. It keeps air in the bottle and away from baby's tummy to help reduce gas, spit-up and burping.

## Less fussing



Philips Avent anti-colic bottle reduces fussing. Babies fed with Philips Avent anti-colic bottles experienced 60% less fussing at night, than babies fed with a leading competitor's vented bottle.\*

## Uninterrupted feeding



The nipple shape allows secure latch and the ribbed texture helps prevents collapse for uninterrupted, comfortable feeding.

## Leak-free

Our anti-colic bottle is designed to prevent leakage whilst feeding, for a truly enjoyable feeding experience.

## Secure latch

Nipple shape designed for secure latch

## Easy to assemble



Our anti-colic bottle has fewer parts for quick and simple assembly.

## Growing up with your baby



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

## Orthodontic Pacifier



Flat, drop-shaped symmetrical nipples respect the natural development of your baby's palate, teeth and gums, even if the pacifier ends up upside-down in the mouth.

## Easy to clean



With its wide neck, our bottle is easier to clean especially with our curved brush head and molded handle-tip that reach the corners, nipples and feeding equipment for thorough cleaning. Ease of cleaning helps to ensure more hygienic feeding.

## BPA free



Philips Avent anti-colic bottle is made of BPA free material (PP).

## Compatible with Classic+ range



Full compatibility with Classic+ bottle range

# Specifications

## Country of origin

**Assembled in UK or in the USA:** With components from the USA, UK and Malaysia

## Material

**Bottle:** Polypropylene, BPA free  
**Nipple:** Silicone, BPA free

## What is included

**Bottle and nipple brush:** 1 pcs  
**Milk powder dispenser:** 1 pcs  
**Translucent pacifier:** 1 pcs  
**Trainer cup:** 1 pcs  
**Baby bottle:** 4 pcs

## Bottle

**Material:** BPA free

## Design

**Bottle design:** Ergonomic shape, Wide neck

## Ease of use

**Bottle use:** Dishwasher & microwave safe  
**Ease of use:** Easy to assemble, Easy to clean, Easy to hold

## Development stages

**Stages:** 0-12 months

## Functions

**Venting system:** Airflex venting system designed to reduce air ingestion

**Nipple:** Easy latch on, Ribbed texture prevents nipple collapse, Proven anti-colic system

**Ease of use:** easy to clean and assemble, leak-free, 4 pieces for easy assembly

**Compatibility:** Full Classic+ compatibility

