

PHILIPS

AVENT

Single Electric Breastfeeding set

Single

Express, store and feed.

Corded use

SCD395/21



Express, Store and Feed Set

Natural motion technology for a quicker milk flow*

The ideal set for mums who plan to express; with all you need to express and store your breast milk and feed your baby. Our baby bottles with an ultra-soft teat closely resemble the breast, making it easy to combine bottle and breast feeding.

Express

- Natural motion technology for a quick milk flow*
- One-size soft and adaptive silicone cushion
- Express without leaning forward
- Personalised experience, 8 + 16 setting levels
- Quiet motor for a discreet experience anywhere (corded use)

Store and transport securely


- For safe storage and transport
- Convenient adapters included
- For use in fridge and freezer

Easy to combine breast and bottle feeding

- Natural latch-on due to the wide breast shaped teat
- Ultra soft teat designed to mimic the feel of the breast
- Anti-colic valve designed to reduce colic and discomfort


Highlights

Quick milk flow*



Express more milk in less time* with a cushion that stimulates the breast to express milk just like baby. It seamlessly adjusts from stimulation mode to expression mode and applies just the right amount of teat stimulation and suction for maximum milk flow. Based on milk flow initiation time (time to Milk Ejection Reflex — MER) results.*

Adapts to your nipple



One size fits all. Because we all come in different shapes and sizes, the silicone cushion gently flexes and adapts to fit your nipple. It fits 99.98% of nipple sizes* (up to 30 mm).

Sit comfortably



Feel relaxed as you express thanks to a design that lets you sit upright rather than having to lean forward. No need to lean forward, a clinically proven comfortable expressing position*. Based on product clinical trial results with 20 participants (2019); 90% of participants find the expressing position comfortable (single electric); 95% of participants find the expressing position comfortable (double electric).

Personalise to your needs



Finely tune every session to your needs with a wide range of stimulation and expression settings. Our breast pump offers 8 stimulation and 16 expression levels for a personalised experience.

Express anywhere discreetly



Express whenever and wherever you like thanks to our quiet motor with a slimline design. Corded use only.

Leak-proof lid



Philips Avent storage cups have a lid that creates a secure seal for safe storage and transport.

Convenient adapters included



For expressing and feeding.

For use in fridge and freezer



For maximum flexibility.

Natural latch on



The wide breast shaped teat promotes a natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Ultra soft teat



The teat has an ultra soft texture, designed to mimic the feel of the breast.

Anti-colic valve



Anti-colic valve designed to keep air away from your baby's tummy, to help reduce colic and discomfort.

Specifications

Power

Mains voltage: 100 V - 240 V

Functions

No leaning forward: Sit in a comfortable position

Soft and adaptive cushion: Gentle stimulation

Settings: 16 expression levels, 8 stimulation levels

Ease of use

Breast pump use: Easy cleaning and assembly, quiet motor, memory function

Material

Bottle: BPA-free*, Polypropylene

Teat: BPA-free*, Silicone

Breast pump: BPA free* (food contact parts only)

What is included

Expression kit: 1 pcs

Motor unit (Corded use only): 1 pcs

120 ml (4 oz) Bottle with Om+ nipple: 1 pcs

9 oz. Bottle with 1 m+ teat: 1 pcs

Micro-USB adapter (corded use): 1 pcs

Milk storage cup (180 ml/6 oz): 5 pcs

Milk storage cup adapter: 1 pcs

Insulation pouch: 1 pcs

Sealing disc: 2 pcs

Disposable breast pads: 2 pcs



* Based on milk flow initiation time (time to Milk Ejection Reflex — MER) results from clinical trial with 20 participants (Netherlands, 2019) compared to time to MER results for other Philips predecessor pump technology from Feasibility study with 9 participants (Netherlands, 2018)

* Based on time to MER results for 1k cushion from clinical trial with 20 participants (Netherlands, 2019) compared to time to MER results from other Philips breast pump in Feasibility study with 9 participants (Netherlands, 2018)

* 1)Mangel et al. Breastfeeding difficulties, breastfeeding duration, maternal body mass index, and breast anatomy: are they related?. Breastfeeding Medicine, 2019, (109 participants, Israel); (2)Ziemer et al. Skin changes and pain in the nipple during the 1st week of lactation.

* Journal of Obstetric, Gynecologic & Neonatal Nursing, 1993, (20 Caucasian participants, USA); (3)Ramsay et al. Anatomy of the lactating human breast redefined with ultrasound imaging, 2005, (28 participants, Australia).

* BPA-free breast pump: Only associated with the bottle and other parts that come into contact with breast milk. Following EU regulation, 10/2011