

PHILIPS

AVENT

Anti-colic bottle gift set

All in One set

SCD397/01



Clinically proven to reduce colic and discomfort*

Designed to reduce gas

The Philips Avent All in One set is ideal for travel, storage, feeding, cleaning and sterilizing baby products.

Other benefits

- This bottle is BPA free
- Full compatibility with Classic+ bottle range
- Compatible range from breastfeeding to cup
- Dispenser holds enough formula for three 260ml/9oz feeds
- Orthodontic, symmetrical and collapsible pacifier nipple

Clinically proven to reduce colic and discomfort*

- Anti-colic system proven to reduce colic*
- 60% less fussing at night*
- Nipple shape designed for secure latch
- Ribbed texture prevents collapse for uninterrupted feeding

Easy to use, clean and assemble

- Easy to assemble with only fewer parts
- Leak-free design
- Easy to clean bottles especially with curved brush head


Highlights

Less colic



Our clinically proven Airflex venting system is designed to reduce gas and colic. As your baby feeds, the valve integrated into the nipple flexes to allow air into the bottle to prevent vacuum build up and vents it towards the back of the bottle. It keeps air in the bottle and away from baby's tummy to help reduce gas, spit-up and burping.

Less fussing



Philips Avent anti-colic bottle reduces fussing. Babies fed with Philips Avent anti-colic bottles experienced 60% less fussing at night, than babies fed with a leading competitor's vented bottle. *

Uninterrupted feeding



The nipple shape allows secure latch and the ribbed texture helps prevent collapse for uninterrupted, comfortable feeding.

Leak-free

Our anti-colic bottle is designed to prevent leakage whilst feeding, for a truly enjoyable feeding experience.

Secure latch

Nipple shape designed for secure latch

Easy to assemble



Our Anti-colic bottle has fewer parts for quick and simple assembly.

Growing up with your baby



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

Orthodontic Pacifier



Flat, drop-shaped symmetrical nipples respect the natural development of your baby's palate, teeth and gums, even if the pacifier ends up upside-down in the mouth.

Easy to clean



With its wide neck, our bottle is easy to clean especially with our curved brush head and molded handle-tip that reach the corners, nipples and feeding equipment for thorough cleaning. Ease of cleaning helps to ensure more hygienic feeding.

BPA free



Philips Avent anti-colic bottle is made of BPA free material (PP).

Compatible with Classic+ range

Full compatibility with Classic+ bottle range

Formula Dispenser



The Philips Avent milk powder dispenser holds 3 pre-measured doses of milk powder, and is ideal for travel.

Specifications

Country of origin

Assembled in UK or in the USA: With components from the USA, UK and Malaysia

Material

Bottle: Polypropylene, BPA free

Nipple: Silicone, BPA free

What is included

Baby bottle: 5 pcs

Bottle and nipple brush: 1 pcs

Milk powder dispenser: 1 pcs

Orthodontic pacifier: 1 pcs

Handles for trainer cup: 1 pcs

Microwave steam sterilizer: 1 pcs

Soft spout: 1 pcs

Bottle

Material: BPA free

Design

Bottle design: Ergonomic shape, Wide neck

Ease of use

Bottle use: Dishwasher & microwave safe

Ease of use: Easy to assemble, Easy to clean, Easy to hold

Development stages

Stages: 0-12 months

Functions

Anti-colic valve: Airflex venting system designed to reduce air ingestion

Nipple: Easy latch on, Ribbed texture prevents nipple collapse, Proven anti-colic system

Ease of use: easy to clean and assemble, leak-free, 4 pieces for easy assembly

Compatibility: Full Classic+ compatibility

