

PHILIPS

AVENIT

Anti-colic bottle gift set

Essentials Set

SCD398/01



Clinically proven to reduce colic and discomfort*

Designed to reduce gas

The Philips Avent Essentials set is ideal for travel, storage, feeding, cleaning and sterilizing baby products.

Clinically proven to reduce colic and discomfort*

- Anti-colic system proven to reduce colic*
- 60% less fussing at night*
- Ribbed texture prevents collapse for uninterrupted feeding
- Nipple shape designed for secure latch

Easy to use, clean and assemble

- Leak-free design
- Easy to assemble with only fewer parts
- Easy to clean bottles especially with curved brush head

Other benefits

- Compatible range from breastfeeding to cup
- Orthodontic, symmetrical and collapsible pacifier nipple
- This bottle is BPA free
- Compatible with Anti-colic range
- Dispenser holds enough formula for three 260ml/9oz feeds


Highlights

Less colic



Our clinically proven Airflex venting system is designed to reduce gas and colic. As your baby feeds, the valve integrated into the nipple flexes to allow air into the bottle to prevent vacuum build up and vents it towards the back of the bottle. It keeps air in the bottle and away from baby's tummy to help reduce gas, spit-up and burping.

Less fussing



Philips Avent anti-colic bottle reduces fussing. Babies fed with Philips Avent anti-colic bottles experienced 60% less fussing at night, than babies fed with a leading competitor's vented bottle.*

Uninterrupted feeding



The nipple shape allows secure latch and the ribbed texture helps prevents collapse for uninterrupted, comfortable feeding.

Leak-free

Our anti-colic bottle is designed to prevent leakage whilst feeding, for a truly enjoyable feeding experience.

Secure latch

Nipple shape designed for secure latch

Easy to assemble



Our Anti-colic bottle has fewer parts for quick and simple assembly.

Growing up with your baby



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

Orthodontic Pacifier



Flat, drop-shaped symmetrical nipples respect the natural development of your baby's palate, teeth and gums, even if the pacifier ends up upside-down in the mouth.

Easy to clean



With its wide neck, our bottle is easy to clean especially with our curved brush head and molded handle-tip that reach the corners,

nipples and feeding equipment for thorough cleaning. Ease of cleaning helps to ensure more hygienic feeding.

BPA free



Philips Avent anti-colic bottle is made of BPA free material (PP).

Compatible Anti-colic range



Full Compatibility with Philips Avent Anti-colic bottle with or without AirFree vent.

Formula Dispenser



The Philips Avent milk powder dispenser holds 3 pre-measured doses of milk powder, and is ideal for travel.

Specifications

Material

Bottle: BPA free, Polypropylene
Nipple: BPA free, Silicone

What is included

Bottle and nipple brush: 1 pcs
Handles for trainer cup: 1 pcs
Milk powder dispenser: 1 pcs
Soft spout: 1 pcs

Microwave steam sterilizer: 1 pcs

Orthodontic pacifier: 1 pcs

Baby bottle: 5 pcs

Design

Bottle design: Ergonomic shape, Wide neck

Ease of use

Bottle use: Dishwasher & microwave safe
Ease of use: Easy to assemble, Easy to clean, Easy to hold

Development stages

Stages: 0-12 months

Functions

Anti-colic valve: Airflex venting system designed to reduce air ingestion
Nipple: Easy latch on, Ribbed texture prevents nipple collapse, Proven anti-colic system
Ease of use: easy to clean and assemble, 4 pieces for easy assembly, leak-free
Compatibility: Full Classic+ compatibility

