

PHILIPS

Analogue baby
monitor

SCD450/05



Two parent units

Both parents can constantly stay in touch with the baby. Two parent units for freedom of movement without losing contact with your baby. With small, mobile and rechargeable parent unit, belt clip and digital signal confirmation.

Complete peace of mind

- Minimum interference from other monitors
- Uses lights to indicate if the baby makes a noise
- Indicates a continuous connection between the two units
- Lights up when batteries need recharging or replacing

Ultimate mobility in and around the house

- Now both parents can monitor baby
- Small, mobile parent unit to freely roam around the house
- Ensures secure reception in and around the house
- Rechargeable, cordless monitoring anywhere
- Belt clip and neck cord for hands-free convenience

Convenience for you and your child

- Adjust volume and sound sensitivity to suit any situation

Soothing words

- Soft glow comforts baby if she wakes

Highlights

16 channels

Minimum interference from other monitors.

Small, mobile parent unit



Allows you to freely roam around the house due to the small size, but still hear when your baby needs you.

Up to 200 m range

Ensures secure reception in and around the house.

Cordless parent unit

Rechargeable, cordless monitoring anywhere.

Soothing night-light



Soft glow comforts baby if she wakes

Specifications

Tuner/Reception/Transmission

Frequency band: 40.675 - 40.690 MHz

Pilot tone: Digital

Maximum range: 200 m

Number of channels: 16

Convenience

Auto out of range warning

Battery charging indication

Battery low indication

Power on indication

Sound-level lights

System check: Battery check light, System

check light, System check alarm

Volume control

Accessories

AC/DC adapter: 3 x 9 V DC 200 mA

Batteries: 3 x AAA rech. batt pack, small parent unit (incl.)

Belt clip: Belt clip on both parent units

Charger: Charging base for parent unit

Docking station: For small parent unit

User manual

Quick start guide

Power

Operating time on battery: 12 hour(s)

