



Purple Fashion Gift Set

- Feeding bottles, pacifiers
- 9oz/260ml
- Slow flow nipple
- Natural nipple shape

SCD693/23

The most natural way to bottle feed

Easy to combine with breastfeeding

The Purple Fashion Gift Set is the ideal gift to give to moms and their precious babies. The gift set comes with two Philips Avent Natural feeding bottles and two pacifiers, to give baby a fashionable feeding solution.

Benefits

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple

A more comfortable and contented feed for baby

- Unique Petals for a soft, flexible nipple without collapse

Vents air away from your baby's tummy

- Advanced anti-colic system with innovative twin valve

Easy to hold and grip in any direction

- Ergonomic shape for maximum comfort

Other benefits

- Orthodontic pacifier
- Compatible with the Philips Avent Natural range
- Simple to use and clean, quick and easy assembly
- This bottle is BPA free*

PHILIPS

Features

Natural latch on



The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Unique comfort petals



Petals inside the nipple increase softness and flexibility without nipple collapse. Your baby will enjoy a more comfortable and contented feed.

Advanced anti-colic valve



Innovative twin valve designed to reduce colic and discomfort by venting air into the bottle and not baby's tummy.

Orthodontic pacifier



The orthodontic, collapsible nipple of the translucent pacifier respects the natural development of baby's palate, teeth and gums. All our pacifiers are made from silicone and are taste and odour-free.

Compatibility across the range



The new Philips Avent Natural bottle is compatible with the Philips Avent range, excluding Classic bottles and cup handles. We advise to use the Natural bottles with Natural feeding nipples only.

Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

Ergonomic shape



Due to the unique shape, the bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

This bottle is BPA free*



The Philips Avent Natural bottle is made of BPA free* material (polypropylene).

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Country of origin

Indonesia Yes

Material

Bottle BPA free*
Polypropylene
Nipple BPA free*
Silicone

What is included

Baby bottle 2 pcs
Pacifier 2 pcs

Bottle

Material BPA free*

Design

Bottle design Ergonomic shape
Wide neck

Ease of use

Bottle use Easy to assemble
Easy to clean
Easy to hold

Development stages

Stages 0 - 6 months

Functions

Anti-colic valve
Latch on
Nipple

Advanced anti-colic system
Easy combine breast and bottle
Natural latch on
Extra soft and flexible teat
Unique comfort petals

Features

Can be sterilized Yes

* 0% BPA, following EU regulation 10/2011



data subject to change
2017, July 26

Version: 3.1.1
UPC:

© 2017 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com