

PHILIPS

AVENT

Natural baby bottle

4 Bottles

9oz/260ml

Slow flow nipple

1m+



SCF013/43



The most natural way to bottle feed

Natural latch on

Our Natural bottle with an ultra soft nipple more closely resembles the breast. The wide breast-shaped nipple with flexible spiral design and comfort petals allows natural latch on and makes it easy to combine breast and bottle feeding.

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple
- Soft nipple designed to mimic the feel of the breast

A comfortable and contented feed for your baby

- Flexible spiral design combined with comfort petals

Vents air away from your baby's tummy

- Unique anti-colic valve technology

Other benefits

- Nipples with different flow rates available
- Ergonomic shape for maximum comfort
- Simple to use and clean, quick and easy assembly
- Compatible with the Philips Avent Natural range
- This bottle is BPA free*

Highlights

Natural latch on



The wide breast shaped nipple promotes a natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Ultra soft nipple



The nipple has an ultra soft texture, designed to mimic the feel of the breast.

Flexible spiral design




Flexible spiral design, combined with our unique comfort petals to create a flexible nipple, allowing for a more natural feed without nipple collapse.

Anti-colic valve technology



Our unique anti-colic valve technology was designed to reduce colic and discomfort, by venting air away from your baby's tummy.

Choosing the right nipple



The Philips Avent Natural range offers different nipple softness and increasing flow rates for every development stage of your baby.

Ergonomic shape



Due to the unique shape, the bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

Compatibility across the range



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

This bottle is BPA free*



The Philips Avent Natural bottle is made of BPA free* material (polypropylene).

Specifications

Country of origin

Indonesia

Material

Bottle: Polypropylene, BPA free*

Nipple: Silicone, BPA free*

What is included

Baby bottle: 4 pcs

Design

Bottle design: Ergonomic shape, Wide neck

Ease of use

Bottle use: Easy to hold, Easy to assemble,
Easy to clean

Development stages

Stages: 0-12 months

Functions

Latch on: Natural latch on, Easy combine breast and bottle

Nipple: Unique comfort petals, Ultra soft and flexible nipple

Anti-colic valve: Advanced anti-colic system

