



Natural baby bottle

- 4 Bottles
- 9oz/260ml
- Slow flow nipple
- 1m+

SCF013/43

The most natural way to bottle feed

Natural latch on

Our Natural bottle with an ultra soft nipple more closely resembles the breast. The wide breast-shaped nipple with flexible spiral design and comfort petals allows natural latch on and makes it easy to combine breast and bottle feeding.

Benefits

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple
- Soft nipple designed to mimic the feel of the breast

A comfortable and contented feed for your baby

- Flexible spiral design combined with comfort petals

Vents air away from your baby's tummy

- Unique anti-colic valve technology

Other benefits

- Nipples with different flow rates available
- Ergonomic shape for maximum comfort
- Simple to use and clean, quick and easy assembly
- Compatible with the Philips Avent Natural range
- This bottle is BPA free*

PHILIPS

Features

Natural latch on



The wide breast shaped nipple promotes a natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Ultra soft nipple



The nipple has an ultra soft texture, designed to mimic the feel of the breast.

Flexible spiral design



Flexible spiral design, combined with our unique comfort petals to create a flexible nipple, allowing for a more natural feed without nipple collapse.

Anti-colic valve technology



Our unique anti-colic valve technology was designed to reduce colic and discomfort, by venting air away from your baby's tummy.

Choosing the right nipple



The Philips Avent Natural range offers different nipple softness and increasing flow rates for every development stage of your baby.

Ergonomic shape



Due to the unique shape, the bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

Compatibility across the range



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

This bottle is BPA free*



The Philips Avent Natural bottle is made of BPA free* material (polypropylene).

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Country of origin

Indonesia Yes

Material

Bottle Polypropylene
BPA free*
Nipple Silicone

data subject to change
2019, July 10

	BPA free*		Easy to clean
What is included		Development stages	
Baby bottle	4 pcs	Stages	0-12 months
Design		Functions	
Bottle design	Ergonomic shape Wide neck	Latch on	Natural latch on Easy combine breast and bottle
Ease of use		Nipple	Unique comfort petals Ultra soft and flexible nipple
Bottle use	Easy to hold Easy to assemble	Anti-colic valve	Advanced anti-colic system

* 0% BPA, following EU regulation 10/2011

* What colic is, and how it affects babies? Colic is caused in part by swallowing air while feeding, which creates discomfort in a baby's digestive system. Symptoms include crying, fussing, gassiness and spit-up.



data subject to change
2019, July 10

Version: 5.0.1
UPC: 00075020068378

© 2019 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com