

PHILIPS

AVENT

Natural baby bottle

3 Bottles

11oz/330ml

Fast flow nipple

6m+

SCF016



The most natural way to bottle feed

Natural latch on

Our bottle with soft and anti-collapse ribbed nipple is designed for growing babies. The comfort petals and natural nipple shape allows natural latch on and makes it easy to combine breast and bottle feeding.

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast shaped nipple
- Soft and smooth silicone for your baby's changing needs

A comfortable and contented feed for your baby

- Flexible, anti-collapse ribbed teat design

Vents air away from your baby's tummy

- Unique anti-colic valve technology

Other benefits

- Nipples with different flow rates available
- Ergonomic shape for maximum comfort
- Simple to use and clean, quick and easy assembly
- Compatible with the Philips Avent Natural range
- This bottle is BPA free*

Highlights

Natural latch on



The wide breast shaped nipple promotes a natural latch on similar to the breast. Designed to give a more comfortable and contented feed for your baby.

Soft and smooth silicone

Bite resistant and smooth nipple designed for your baby's changing needs.

Flexible, anti-collapse design




Petals and ridges inside the nipple increase softness and flexibility without nipple collapse.

Anti-colic valve technology



Our unique anti-colic valve technology was designed to reduce colic and discomfort, by venting air away from your baby's tummy.

Choosing the right nipple



The Philips Avent Natural range offers different nipple softness and increasing flow rates for every development stage of your baby.

Ergonomic shape



Due to the unique shape, the bottle is easy to hold and grip in any direction for maximum comfort, even for baby's tiny hands.

Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

Compatibility across the range



Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it!

This bottle is BPA free*



The Philips Avent Natural bottle is made of BPA free* material (polypropylene).

Specifications

Material

Bottle: Polypropylene, BPA free*

Nipple: Silicone, BPA free*

What is included

Baby bottle: 3 pcs

Bottle

Material: BPA free*

Design

Bottle design: Ergonomic shape, Wide neck

Ease of use

Bottle use: Easy to hold, Easy to assemble,
Easy to clean

Development stages

Stage: 0-12 months

Functions

Latch on: Natural latch on, Easy combine breast
and bottle

Nipple: Unique comfort petals, Soft and anti-
collapse

Anti-colic valve: Advanced anti-colic system

